

PC-EHSP-9-13-21 ©ECEP HAN CATALOGUE LANGUAGE: EN





CREATIVE



HEALTY



TRUSTWORTHY







Ecep Han Safety Products (EHSP $^{TM}$ ) is the brand under which we produce medical, hygienic, medical devices, surgical drapes and Personal Protective Equipment(PPE)

# **ABOUT US**

# Multinational Company Specialized in Medical Devices, Covid-19 PPE and Surgical Drapes Manufacturing

Medical professionals, textile engineers, development experts, and good buisness people gathered from corners of the world to make a single soul and contribute to humanity's health and good. We wake up thinking about you and hit the hay, contemplating what boon we have delivered to you. In world where health is only buisness, in ECEP HAN worldwide, we are putting hands to hands innovating in speed and accuracy to produce medical textiles and PPE where it will go for good.





# Table of Contents

EHSP-MASK	3
EHSP-COVERALLS	19
EHSP-GOWNS	27
EHSP-SCRUBS	37
EHSP-ACCESSORIES	43
EHSP-HYGIENE PRODUC	ΓS59
EHSP-GOWNS GUIDE	65
EHSP-MASK GUIDE	72
EHSP-CERTIFICATES AND	TEST REPORTS81



Medical professionals, textile engineers, development experts, good business and people gathered from corners of the world to make a single soul and contribute to humanity's health and good. We wake up thinking about you and hit the hay, contemplating what boon we have delivered to you. In a world where health is only business, we ECEP HAN continually at innovate medical products with speed and accuracy to deliver health support around the world for good.



Trade
Hear
Speak
Work
&
Communicate

# for Good




# MEDICAL, RESPIRATOR AND CIVILIAN FACE MASKS



# EHSP 300 N95/FFP2 NR RESPIRATOR **FACE MASK**

N95/FFP2 NR Medical Grade Respirator Face Mask





5 Layers Face Mask

Spunbond Meltblown Meltblown Hot cotton Spunbond



Respiratory Protective Devices - Filtering to Protect Against Particles



BFE > 95% PFE < 99.97%



Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059





#### **GENERAL DESCRIPTION**

- N95/FFP2 High-Efficiency Filtration > 95% < 99.97%
- 5 Layers of Protection
- EN 149:2001+A1:2010 Tested
- **NIOSH Tested**
- Elasticated Earloop
- Ultrasonic Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- EN13274-7:2019 Tested at 95 l/min Percentage of Flow Rate for Penetration of Paraffin Oil Mist



Adjustable Nose Bridge

Sterility

upon Request

**Fastening Type** 

Elasticated Earloop

Standard Adult Size





Gender Unisex

Disposable



**Private Labeling** Possible upon Request; Subject to MOQ.



SHIPMENT DETAILS AND SERVICES

**Packing Style** 1 PCs, 50 PCs, 60 PCs

- Health Care Providers
- First Responders
- To Protect Against Airborne Particles, or Liquid Contaminating

## **RECOMMENDED TO:**

- Face





Airway Shipment



Seaway Shipment



Railway Shipment















CFR 2020



(1000 x 1200 mm) (40" x 48") Europe, North America

Standard S



Àsia (800 x 1200 mm)

1100 x 1100 mm)

 $(44" \times 44")$ 



(32" x 48") Europe





Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkev@ecephan.com



ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA








SWISS Office:
ECEP HAN MEDIKAL TEK, SAN, TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND



+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



# EHSP 310 N95/FFP2 NR RESPIRATOR **FACE MASK HEADLOOP STYLE**

N95/FFP2 NR Medical Graded Face Mask





5 Layers Face Mask Pro-Muscle Non-Woven Fabric Meltblown Meltblown Filter Sponge

Non-Woven Fabric



Respiratory Protective Devices - Filtering to Protect Against Particles



BFE> 95% PFE< 99.97%



Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059





#### **GENERAL DESCRIPTION**

- N95/FFP2 High-Efficiency Filtration > 95% < 99.97%
- 5 Layers of Protection
- EN 149:2001 + A1:2010 and NIOSH Tested
- Headloop Style
- Elasticated Headloops
- Ultrasonic Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- E13274-7:2019 Tested at Flow Rate Percentage of 95 I/min for Penetration of Paraffin Oil Mist

Adjustable Nose Bridge



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Elasticated Headloops



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Adult Size



SHIPMENT DETAILS AND SERVICES

**Packing Style** 1 PC, 50 PCs, 60 PCs

- Health Care Providers
- First Responders

Innovation.Speed.Accuracy

To Protect Against Airborne

- Particles, or Liquid Contaminating Face

#### **RECOMMENDED TO:**



Roadway Shipment

Airway Shipment

Seaway Shipment



Railway Shipment

#### 🥌 🚤 🚜 ž INCOTERMS









CFR 2020



(40" x 48") Europe, North America

Standard S



1100 x 1100 mm)  $(44" \times 44")$ Àsia



(800 x 1200 mm) (32" x 48") Europe



TURKEY O'MCE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA






SWISS Office:
ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND



+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY





# EHSP 320 N95/FFP2 NR RESPIRATOR **FACE MASK**

N95/FFP2 NR Medical Grade Respirator Face Mask





5 Layers Face Mask

Spunbond Meltblown Meltblown Hot cotton Spunbond



Respiratory Protective Devices - Filtering to Protect Against Particles



BFE > 95% PFE < 99.97%



Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059





#### **GENERAL DESCRIPTION**

- N95/FFP2 High-Efficiency Filtration > 95% < 99.97%
- 5 Layers of Protection
- EN 149:2001+A1:2010 Tested
- **NIOSH Tested**
- Elasticated Earloop
- Ultrasonic Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- EN13274-7:2019 Tested at 95 l/min Percentage of Flow Rate for Penetration of Paraffin Oil Mist



Adjustable Nose Bridge



Disposable







Gender Unisex



**Fastening Type** Elasticated Earloop



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Adult Size

+ 1 +

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



SHIPMENT DETAILS AND SERVICES

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

**Packing Style** 1 PCs, 50 PCs, 60 PCs

- Health Care Providers
- First Responders

Innovation.Speed.Accuracy

 To Protect Against Airborne Particles, or Liquid Contaminating

#### **RECOMMENDED TO:**

- Face



ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





SWISS Office:
ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Standard S

(1000 x 1200 mm)

1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Europe

Àsia

Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



dubai@ecephan.com kabul@ecephan.com

# EHSP 330 N95/FFP2 NR RESPIRATOR **FACE MASK**

N95/FFP2 NR Medical Grade Respirator Face Mask





5 Layers Face Mask

Spunbond Meltblown Meltblown Hot cotton Spunbond



Respiratory Protective Devices - Filtering to Protect Against Particles



BFE > 95% PFE < 99.97%



Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059





#### **GENERAL DESCRIPTION**

- N95/FFP2 High-Efficiency Filtration > 95% < 99.97%
- 5 Layers of Protection
- EN 149:2001+A1:2010 Tested
- **NIOSH Tested**
- Elasticated Earloop
- Ultrasonic Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- EN13274-7:2019 Tested at 95 l/min Percentage of Flow Rate for Penetration of Paraffin Oil Mist



Adjustable Nose Bridge

Non-Sterile; Sterilized

upon Request



Disposable



Gender Unisex



**Private Labeling** Possible upon Request; Subject to MOQ.



SHIPMENT DETAILS AND SERVICES

**Packing Style** 1 PCs, 50 PCs, 60 PCs

#### **RECOMMENDED TO:**

- Health Care Providers
- First Responders
- To Protect Against Airborne Particles, or Liquid Contaminating Face



**Fastening Type** Elasticated Earloop

Sterility



Standard Adult Size

## 🥌 🚤 🚜 ž + 1 +



Seaway Shipment



Roadway Shipment



Railway Shipment















(1000 x 1200 mm) (40" x 48") Europe, North America

Standard S



Àsia (800 x 1200 mm)

 $(44" \times 44")$ 

1100 x 1100 mm)



(32" x 48") Europe









ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





SWISS Office:
ECEP HAN MEDIKAL TEK, SAN, TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND



swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





+93 202 501 610

# EHSP 340 N95/FFP2 REUSABLE AND **WASHABLE FACE MASK**

N95/FFP2 Reusable, washable high-quality face mask





3 Layers Fabric Obtained from Anti-Microbial Yarn by Adding Silver Ions Into the Fiber Polymer Structure



Re-Usable Washable Durable



BEE> 95%





#### **GENERAL DESCRIPTION**

- N95/FFP2 High-Efficiency Filtration > 95% < 99.97%
- COVID-19 Protection
- Filtration has been Tested
- Elasticated Earloop
- Soft and Comfortable Edges
- Individually Tailored
- Color Options
- Eco Friendly
- Durable
- **UV** Disinfected
- Small, Large and Medium for Better Fit

**Microbial Cleanliness** EN 14683+AC:2019



Reusable



Sterility Non-Sterile; Sterilized upon Request



**UV** Disinfected Each Mask is UV Disinfected



**Equipped with FFP2 Valve** Protects Solid Dust & Oil Mist (EN 149:2001 + A1:2009)



**Private Labeling** Possible upon Request; Subject to MOQ.



Different Adult Size Small, Medium, Large



**Packing Style** 10, 50 PCs, and Individually Packed

#### **RECOMMENDED TO:**

- Normal Outings
- First Responders
- To Protect Against Airborne Particles, or Liquid Contaminating Face

Airway Shipment



Seaway Shipment



Roadway Shipment



Railway Shipment







SHIPMENT DETAILS AND SERVICES







EXW 2020 FOB 2020



CIF 2020 DAP 2020



CFR 2020



(1000 x 1200 mm) (40" x 48") Europe, North America

Standard



 $(44" \times 44")$ Àsia

1100 x 1100 mm)



(800 x 1200 mm) (32" x 48") Europe





Turkey Office: ECEP HAN MEDİKAL TEK, SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA



+1 480 656 5519













+93 202 501 610

# **EHSP 410 DISPOSABLE S-MASK**

Disposable Level-1 (Type-I) S-Mask with Option of Normal Mask Delta P < 5.0 mm H<sub>2</sub>O/cm<sup>2</sup>





3 Ply S-Mask Spunbond Meltblown Spunbond



Sub-Micron Particulate Filtration



(BFE > 95%, PFE > 95% @ 0.1micron)



Fluid Penetration Resistance (80 mmHg)





#### **GENERAL DESCRIPTION**

- 3 PLY Rectangular S-Mask
- Elasticated Earloop
- Elasticated Neckband 600 mm Long
- Ultrasonic Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- ASTM F2101 (Sub-Micron Particulate Filtration Efficiency) Tested
- ASTM F2299 (BFE> 95%, PFE> 95% @ 0.1micron) Tested
- ASTM F1862 (80 mmHg)
- Flame Spread Class-I

Adjustable Nose Bridge



Disposable



Sterility Non-Sterile; Sterilized upon Request (UV- Sterilized Sterilized)



Gender



Fastening Type Elasticated Earloop & Elasticated Neckband



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Adult Size

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



SHIPMENT DETAILS AND SERVICES

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

Packing Style 10 PCs, 25 PCs & 50 PCs Boxes & Plastic Bags

- Low-Risk Clinical Applications that do not
- Enhancing Infection Control
- Civilian Face-piece



- Preventing Risk of Cross-Contamination

Innovation.Speed.Accuracy

#### **RECOMMENDED TO:**



Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA




SWISS OFFICE:
ECEP HAN MEDIKAL TEK, SAN, TIC LTD, STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Standard

(1000 x 1200 mm)

1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Europe

Àsia





# EHSP 415 DISPOSABLE S-MASK (CHILD)

Disposable Level-II S-Mask, Not Sterilized, Delta P < 5.00 mm H<sub>2</sub>O/cm<sup>2</sup>





3 Ply S-Mask Spunbond Meltblown Spunbond



Sub-Micron Particulate Filtration



(BFE> 98%, PFE> 98% @ 0.1micron)



Fluid Penetration Resistance (120 mmHg)





#### **GENERAL DESCRIPTION**

- PLY Rectangular S-Mask
- Elasticated Earloop
- Elasticated Neckband 400 mm Long
- Ultrasonic Banded Edges
- Seamless Crimping
- Optional Adjustable Nose Bridge
- ASTM F2100 (Sub-Micron Particulate Filtration Efficiency) Tested
- ASTM F2299 (BFE ≥ 98%, PFE ≥ 98% @ 0.1micron) Tested
- Fluid Penetration Resistance (120 mmHg)
- Flame Spread Class-I

Adjustable Nose Bridge upon Request



**Disposable** 



Sterility Non-Sterile; Sterilized upon Request (UV- Sterilized Sterilized)



Gender



Fastening Type Elasticated Earloop & Elasticated Neckband



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard 4-14 Year Children

+ 11 + 111

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



SHIPMENT DETAILS AND SERVICES

🥌 🛶 🚁 💈

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

**Packing Style** 10 PCs, 25 PCs & 50 PCs Plastic Bags & Boxes

- Low-Risk Clinical Applications that do not
- Enhancing Infection Control
- Preventing Risk of Cross-Contamination
- Civilian Face-piece

Innovation.Speed.Accuracu



## **RECOMMENDED TO:**



Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA






SWISS OFFICE:
ECEP HAN MEDIKAL TEK, SAN, TIC LTD, STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354



Standard S

(1000 x 1200 mm)

1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Europe

Àsia

Afghanistan: ECEP HAN TRADING LLC





10

# **EHSP 420 DISPOSABLE S-MASK**

Disposable Level-2 (Type-II) S-Mask with Option of Normal Mask Delta P < 6.0mm H<sub>2</sub>O/cm<sup>2</sup>





3Ply S-Mask Spunbond Meltblown Spunbond



Sub-Micron Particulate Filtration



(BFE > 98%, PFE > 98% @ 0.1micron)



Fluid Penetration Resistance (120 mmHg)





#### **GENERAL DESCRIPTION**

- 3PLY Rectangular S-Mask
- Elasticated Earloop
- Elasticated Neckband 600 mm Long
- Ultrasonic Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- ASTM F2100 (Sub-Micron Particulate Filtration Efficiency) Tested
- ASTM F2299 (BFE> 98%, PFE> 98% @ 0.1micron) Tested
- Fluid Penetration Resistance (120 mmHg)
- Flame Spread Class-I

Adjustable Nose Bridge



Disposable



Non-Sterile; Sterilized upon Request



Gender



Fastening Type Elasticated Earloop & Elasticated Neckband



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Adults Size 6.9" x 3.7" (17.3cm x 9.5cm)

+ 11 + 111



SHIPMENT DETAILS AND SERVICES

**Packing Style** 10 PCs, 25 PCs & 50 PCs Boxes & Plastic Bag

## **RECOMMENDED TO:**

- Low-Risk Clinical Applications that do not Involve Blood-Borne Pathogens or Bodily Fluids
- Enhancing Infection Control
- Preventing Risk of Cross-Contamination
- Civilian Facepiece



Airway Shipment



Railway Shipment

# 🥌 🛶 🚁 💈











CFR 2020



(40" x 48") Europe, North America

1100 x 1100 mm)



 $(44" \times 44")$ Àsia



(800 x 1200 mm) (32" x 48") Europe



+90 212 852 3447 turkev@ecephan.com



ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com



SWISS OFFICE:
ECEP HAN MEDIKAL TEK, SAN, TIC LTD, STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND

+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354

+97 150 696 1828 \ +97 155 676 1828 dubai@ecephan.com



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

+93 202 501 610

# **EHSP 430 DISPOSABLE S-MASK**

Disposable Level-3 (Type-IIR) S-Mask with Option of Normal Mask Delta P < 6.0 mm H<sub>2</sub>O/cm<sup>2</sup>





3 Ply S-Mask Spunbond Meltblown Spunbond



Sub-Micron Particulate Filtration



(BFE> 98%, PFE> 98% @ 0.1micron)



Fluid Penetration Resistance (160 mmHg)





#### **GENERAL DESCRIPTION**

- 3 PLY Rectangular S-Mask
- Elasticated Earloop
- Elasticated Neckband 600 mm Long
- Ultrasonic Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- ASTM F2100 (Sub-Micron Particulate Filtration Efficiency) Tested
- ASTM F2299 (BFE> 98%, PFE> 98% @ 0.1micron) Tested
- Fluid Penetration Resistance (160 mmHg)
- Flame Spread Class-I

Adjustable Nose Bridge



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Elasticated Earloop & Elasticated Neckband



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Adult Size

+ 🚄 + 🛍 + 📖

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



SHIPMENT DETAILS AND SERVICES

G - 2 2 2

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

Packing Style 10 PCs, 25 PCs & 50 PCs Plastic Bag

#### **RECOMMENDED TO:**

Enhancing Infection Control

Innovation.Speed.Accuracy

- Preventing Risk of Cross-Contamination
- Procedures with Moderate to Heavy Amounts of Blood, Fluid, Spray and/or

- Aerosol Exposure



Turkey Office: ECEP HAN MEDİKAL TEK, SAN, TİC LTD. ŞTİ ÇAKMAKLI MAH, DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Standard -

(1000 x 1200 mm)

1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Europe

Àsia

Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



12

# **EHSP 440 SURGICAL TIE-ON MASK**

Disposable Level-2 or 3 (Type-IIR) S-Mask with Option of Normal Mask Delta P < 6.0 mm H<sub>2</sub>O/cm<sup>2</sup>





3 Ply Tie-On Mask 4 Ply upon Request Spunbond

Meltblown Spunbond



Respiratory Protective Devices - Filtering to Protect Against Particles



BFE> 98% PFE< 98%



Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059





#### **GENERAL DESCRIPTION**

- 3 PLY SURGICAL TIE-ON,
- Ultrasonic Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- ASTM Efficiency ASTM F2101 and Sub-Micron Particulate Filtration Efficiency
- ASTM F2299 (PFE ≥ 98% @ 0.1 micron)
- ASTM F2101, EN 14683 BFE ≥ 98%
- Fluid Penetration Resistance ASTM F1862 (160 mmHg)
- Flame Spread Class-I (≥ 3.5 seconds) / EN Medical Directive (2007/47/EC, MDD 93/42/EEC)

Adjustable Nose Bridge



Disposable



Sterility

Non-Sterile; Sterilized upon Request (UV- Sterilized Sterilized)



Gender Unisex



**Fastening Type** Tie-on Head Band



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Adult Size



Packing Style 1 PC, 10 PCs, 25 PCs Sterilized Packing

#### **RECOMMENDED TO:**

- Enhancing Infection Control
- Preventing Risk of Cross-Contamination
- Procedures with Moderate to Heavy Amounts of Blood, Fluid, Spray and/or Aerosol Exposure

# Airway Shipment



Seaway Shipment



Roadway Shipment



Railway Shipment

# + 11 +

# G - 2 2 2

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020



SHIPMENT DETAILS AND SERVICES







(40" x 48") Europe, North America



1100 x 1100 mm)  $(44" \times 44")$ Àsia



(800 x 1200 mm) (32" x 48") Europe





ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND



+41 79 714 3173



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



# **EHSP 450 MASK WITH EYE SHIELD**

Disposable Level-3 (Type-IIR) S-Mask with Option of Normal Mask with Anti-Fog Eye-Shield, Delta P < 6.0 mm H<sub>2</sub>O/cm<sup>2</sup>





3 Ply Tie-On Mask

Spunbond

Meltblown

Spunbond

Anti-Fog PC Protective Film



Sub-Micron Particulate Filtration



(BFE> 98%, PFE> 98% @ 0.1micron)



Fluid Penetration Resistance (160 mmHg)





#### **GENERAL DESCRIPTION**

- 3 PLY Level-3 with Protective Shield
- Banded Edges
- Seamless Crimping
- Adjustable Nose Bridge
- ASTM Efficiency ASTM F2101 and Sub-Micron Particulate Filtration Efficiency
- ASTM F2299 (PFE ≥ 98% @ 0.1micron)
- ASTM F2101, EN 14683 BFE ≥ 98%
- Fluid Penetration Resistance ASTM F1862 (160 mmHg)
- Flame Spread Class-I (≥ 3.5 seconds) / EN Medical Directive (2007/47/EC, MDD 93/42/EEC)

Adjustable Nose Bridge



Sterility

Non-Sterile; Sterilized upon Request (UV- Sterilized Sterilized)



**Fastening Type** Tie-on Head Band



Standard Adult Size

+ 11 +

Airway Shipment

Seaway Shipment

Roadway Shipment



Disposable



Gender Unisex



Private Labeling Possible upon Request; Subject to MOQ.



SHIPMENT DETAILS AND SERVICES

G - 2 2 2

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

Packing Style 1 PC, 10 PCs, 25 PCs Sterilized Packing

#### **RECOMMENDED TO:**

- Enhancing Infection Control
- Preventing Risk of Cross-Contamination
- Procedures with Moderate to Heavy Amounts of Blood, Fluid, Spray and/or

- Aerosol Exposure





USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Standard S

(1000 x 1200 mm)

1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Europe

Àsia

Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY





+90 212 852 3447 turkev@ecephan.com

Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

14

# **EHSP PPP05 FACE SHIELD FOR KIDS**

Anti-Fog Visor, Gripper Mechanism, Polycarbonate Plus 1 Spare Visors





Outer Layer:

100% Original PET Resin

Middle Layer:

100% Original PET Resin

or Recycled PET



Protects Against Splash, Spray and Fluids





#### **GENERAL DESCRIPTION**

- EHSP-PPP05 Face Shields are Designed to Prevent the Splash, Spray and Transfer of Non-Corrosive Fluids
- Adjustable Rubber Band to Fit Head Size
- Adhesive Tape on the Sponge to Protect Forehead
- Crown Exterior Size: 17.3"\*0.98" (46\*2.5 cm)
- Crown Interior Size: 15.5"\*0.59" (38.5\*1.5 cm)
- Visor Size: 12.6"\*7.48" (32\*19 cm)
- Elongation: 4%
- Tensile Strength: 47 MPa
- Impact Resistance: 595 KJm<sup>2</sup>
- Density: 1.35 g/cm<sup>3</sup> ± 0.05 g/cm<sup>3</sup>
- Melting Point: 243°C
- Surface Tension: ≥38 dyn

Visor Thickness: 150µ-800µ Cut Dimensions: ≤1200\*1600 mm



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Elasticated Headloop, Adjustable Rubber Band



**Private Labeling** Possible upon Request; Subject to MOQ.



Size Child Size



**Packing Style** 50 PCs Individually Packed With an Extra Visor

## **RECOMMENDED TO:**

- Health Worker
- People Getting Very Close with Others Who Can be Infected
- Be Used with a Mask for Optimum Results

#### SHIPMENT DETAILS AND SERVICES

#### G - 2 2 2 + 11 + INCOTERMS







Roadway Shipment



Railway Shipment











CFR 2020



(40" x 48") Europe, North America

Standard -



1100 x 1100 mm)  $(44" \times 44")$ Àsia



(800 x 1200 mm) (32" x 48")



Europe



C\*

Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkev@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com

+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354

+97 150 696 1828 \ +97 155 676 1828 dubai@ecephan.com




Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

+93 202 501 610

# **EHSP PPP05 FACE SHIELD**

Anti-Fog Visor, Gripper Mechanism, Polycarbonate Plus 1 Spare Visors





Outer Layer:

100% Original PET Resin

Middle Layer:

100% Original PET Resin

or Recycled PET



Protects Against Splash, Spray and Fluids





#### **GENERAL DESCRIPTION**

- EHSP-PPP05 Face Shields are Designed to Prevent the Splash, Spray and Transfer of Non-Corrosive Fluids
- Adjustable Rubber Band to Fit Head Size
- Adhesive Tape on the Sponge to Protect Forehead
- Crown Exterior Size:17.3"\*0.98" (46\*2.5 cm)
- Crown Interior Size: 15.5"\*0.59" (38.5\*1.5 cm)
- Visor Size: 12.6"\*7.87" (32\*20 cm)
- Elongation: 4%
- Tensile Strength: 47 MPa
- Impact Resistance: 595 KJm<sup>2</sup>
- Density: 1.35 g/cm<sup>3</sup> ± 0.05 g/cm<sup>3</sup>
- Melting Point: 243°C
- Surface Tension: ≥38 dyn



Visor Thickness: 150µ-800µ Cut Dimensions: ≤1200\*1600 mm



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Elasticated Headloop, Adjustable Rubber Band



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Adult Size



**Packing Style** 50 PCs Individually Packed With an Extra Visor

#### **RECOMMENDED TO:**

- Health Worker
- People Getting Very Close with Others Who Can be Infected
- Be Used with a Mask for Optimum Results

# SHIPMENT DETAILS AND SERVICES







Seaway Shipment

Roadway Shipment

Railway Shipment















CFR 2020



(40" x 48") Europe, North America 1100 x 1100 mm)

Standard S



 $(44" \times 44")$ Àsia



(800 x 1200 mm) (32" x 48") Europe





+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND +41 79 714 3173

swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



# **EHSP PROTECTIVE FACE SHIELD WITH GOGGLE**

Ecep Han
Trading for good.

Anti-Fog Visor, Goggle , Polycarbonate Plus 1 Spare Visors



Outer Layer:

100% Original PET Resin

Middle Layer:

100% Original PET Resin

or Recycled PET



Protects Against Splash, Spray and Fluids





#### **GENERAL DESCRIPTION**

 EHSP-Face shield with google are Designed to Prevent the Splash, Spray and Transfer of Non-Corrosive Fluids

Visor Size: 12.6"\*7.48" (32\*19 cm)

Elongation: 4%

Tensile Strength: 47 MPa

Impact Resistance: 595 KJm<sup>2</sup>

Density: 1.35 g/cm<sup>3</sup> ± 0.05 g/cm<sup>3</sup>

Melting Point: 243°C

Surface Tension: ≥38 dyn

Colour: White / Transparent



# Visor Thickness: 150µ-800µ

Cut Dimensions: ≤1200\*1600 mm



**Disposable** 



#### Sterility

Non-Sterile; Sterilized upon Request



Gender Unisex



#### **Fastening Type**

Vented foam design for increased air flow



Private Labeling Possible upon Request; Subject to MOQ.



Size Adult and Child



**Packing Style** Face protection goggle / box

Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to the end-users to make their own determination as to the suitability of the product for their own intended purpose(s)

#### Storage & Handling

Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture and excessive heat above 100°F (37°C).

#### **RECOMMENDED TO:**

- Health Worker
- People Getting Very Close with Others Who Can be Infected
- Be Used with a Mask for Optimum Results

# SHIPMENT DETAILS AND SERVICES

# + 11 +























(1000 x 1200 mm) (40" x 48") Europe, North America

Standard S



1100 x 1100 mm)  $(44" \times 44")$ Àsia



(800 x 1200 mm) (32" x 48")



Europe







+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com

+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354

+97 150 696 1828 \ +97 155 676 1828 dubai@ecephan.com



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



# Protect

Your Physical Health

by getting our good products

Your Mental
Health

by working with our good team.



# **EHSP 600 PROTECTIVE CAT-III TYPE 5B-6B**

AAMI Level-2 CAT-III Type 5B/6B Coverall





(PP+PE) 55-60 **GSM** Laminated



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms



AAMI Level-II Tested





### **GENERAL DESCRIPTION**

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases, Micro-organisms such as Fluids, Aerosols, or Solid Dust Particles
- Not Stipulated by Seam Technology
- Tight Seam Sewn; All Seams Overlocked, Elasticated on Face, Wrist, Rear, and Ankles
- The Zipper is Furnished with a Metal Stopper and Flap Stickers



Breathable, Lightweight, Flexible, Particle-Tight & Spray-Tight Material



**Disposable** 



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Front Zipper with Lap Sticker, Elasticated over Face, Wrists & Waistband



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



SHIPMENT DETAILS AND SERVICES

4 - a a

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

**Packing Style** 32 PCs & 65 PCs

- Healthcare Professionals. Pharmaceutical Industries
- Emergency Wards, Laboratories,
- Outpatients

#### **RECOMMENDED TO:**

- and Isolated Medical Wards

- First Responders, Ambulances Team
- Inspection Room, & other Public Places





C\*

CECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH, DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA









Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUF 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Standard

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Europe

Àsia

(1000 x 1200 mm)

1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America







(PP+PE) 55-60 **GSM** Laminated



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms



AAMI Level-II Tested





#### **GENERAL DESCRIPTION**

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases. Micro-organisms such as Fluids, Aerosols, or Solid Dust Particles
- Not Stipulated by Seam Technology
- Tight Seam Sewn; All Seams Overlocked, Elasticated on Face, Wrist, Rear, and Ankles
- The Zipper is Furnished with a Metal Stopper and Flap Sticker
- +EHSP 500 Built-in Non-Slip Sole Over-locked



Breathable, Lightweight, Flexible, Particle-Tight & Spray-Tight Material



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Front Zipper with Lap Sticker, Elasticated Over Face, Wrists & Waistband



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



SHIPMENT DETAILS AND SERVICES

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

Packing Style 25 PCs & 50 PCs

- Healthcare Professional's Critical Airborne Zones, Limited Hazardous Chemical Splash and Dust Zones
- Protection Against Asbestos, Acid, and Alkali Splash and Contaminated Water
- Industrial Protection Usage

#### **RECOMMENDED TO:**





Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKCEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101 SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUF 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Standard

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Europe

Àsia

(1000 x 1200 mm)

1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America





# EHSP 610 REINFORCED CAT-III TYP 3/4

AAMI Level-3 CAT III Type 3/4 Reinforced Coverall





(PP+PE)



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms



AAMI Level-3 Tested



55-60 **GSM** 



### **GENERAL DESCRIPTION**

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases, Micro-organisms such as Fluids, Aerosols, or Solid Dust Particles
- Not Stipulated by Seam Technology
- Tight Seam Sewn; all Seams Overlocked, Elasticated on Face, Wrist, Rear, and Ankles
- The Zipper is Furnished with a Metal Stopper and Flap Sticker



Breathable, Lightweight, Flexible, Particle-Tight & Spray-Tight Material



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Front Zipper with Lap Sticker, Elasticated over Face, Wrists & Waistband



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL

+ 1 +



SHIPMENT DETAIL

**Packing Style** 32 PCs & 65 PCs

- Pharmaceutical Processing Centers
- Handling Toxic Powders with Particles above 1 micron

- Healthcare Professional's Critical Airborne Zones
- Industrial Protection Usage

# **RECOMMENDED TO:**

# Seaway Shipment

Roadway Shipment

Airway Shipment



Railway Shipment











CFR 2020



(40" x 48") Europe, North America



1100 x 1100 mm)  $(44" \times 44")$ Àsia

Standard



(800 x 1200 mm) (32" x 48") Europe





Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKCEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101 SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUF 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354

+97 150 696 1828 \ +97 155 676 1828 dubai@ecephan.com



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



# **EHSP 610B REINFORCED CAT-III TYPE 3/4**

AAMI Level-3 CAT-III Type 3/4 Reinforced Coverall + EHSP 500 Built-in Boot-Cover (Bunny Style)





(PP+PE) 55-60

**GSM** Laminated



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms



AAMI Level-III Tested





#### **GENERAL DESCRIPTION**

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases. Micro-organisms such as Fluids, Aerosols, or Solid Dust Particles
- Not Stipulated by Seam Technology
- Tight Seam Sewn; All Seams Overlocked, Elasticated on Face, Wrist, Rear, and Ankles
- The Zipper is Furnished with a Metal Stopper and Flap Sticker
- +EHSP 500 Built-in Non-Slip Sole **Boot-Cover**



Breathable, Lightweight, Flexible, Particle-Tight & Spray-Tight Material



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Front Zipper with Lap Sticker, Elasticated over Face, Wrists & Waistband



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL



**Packing Style** 40 PCs & 80 PCs

#### **RECOMMENDED TO:**

- Healthcare Professional's Critical
- Disaster Red-Zones with Particles above 1 micron and Handling Organic and Inorganic Chemicals
- Biological Field Suitable

- Airborne Zones





Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment
































SHIPMENT DETAILS AND SERVICES







Standard



1100 x 1100 mm)  $(44" \times 44")$ Àsia



(800 x 1200 mm) (32" x 48") Europe





Turkey Office: TGIREY OHRE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101 SCOTTSDALE AZ 85260, ARIZONA USA













Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



# **EHSP 650-IND INDUSTRIAL COVERALL**

Disposable Industrial Coverall





40 GSM MD(32) CD(20)



Protects Against Specific Industrial Hazards





#### **GENERAL DESCRIPTION**

- European Standards Certified to Help Protect Against Specific Industrial Hazards
- Comfortable and Easily Wearable
- Recommended for Commercial Buildings, Construction Sites, Facility Sanitation, General Manufacturing, and Food Processing Facilities



The Hood Free Coverall is White in Color with Elastic Cuff, Waist, & Ankle



Disposable



Sterility Non-Sterile



Gender Unisex



Fastening Type Two-Way Zipper, Overlocked Inside, Elasticated Cuff, Waist, and Ankle



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL



Packing Style 40 PCs & 80 PCs

#### **RECOMMENDED TO:**

- Protection Against Chemical Compounds, Dust and Fibers
- Outdoor Painting Industries
- Varnishing Outdoor Industries
- Low Danger Site Inspection

# SHIPMENT DETAILS AND SERVICES

# 

Seaway Shipment

Roadway Shipment

Railway Shipment















(1000 x 1200 mm) (40" x 48") Europe, North America



1100 x 1100 mm) (44" x 44") Àsia

Standard



(800 x 1200 mm) (32" x 48") Europe





Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101 SCOTTSDALE AZ 85260, ARIZONA USA











Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354

+97 150 696 1828 \+97 155 676 1828 dubai@ecephan.com



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



# **EHSP 660-IND INDUSTRIAL COVERALL**

Disposable Industrial Coverall





40 GSM MD(32) CD(20)



Protects Against Specific Industrial Hazards





#### **GENERAL DESCRIPTION**

- European Standards Certified to Help Protect Against Specific Industrial Hazards
- Comfortable and Easily Wearable
- Recommended for Commercial Buildings, Construction Sites, Facility Sanitation, General Manufacturing, and Food Processing Facilities



The Hooded Coverall is White in Color with Elastic Cuff, Waist, & Ankle



Disposable



Sterility Non-Sterile



Gender Unisex



Fastening Type Two-Way Zipper, Overlocked Inside, Elasticated Cuff, Waist, and Ankle



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL



Packing Style 40 PCs & 80 PCs

- Protection Against Chemical Used in Agricultural Processes
- Food Industry
- Low Danger Site Inspection

- Protection Against Chemicals Compounds, Asbestos, Dust & Fibers

#### **RECOMMENDED TO:**

# Airway Shipment Seaway Shipment DAP 2020 Roadway Shipment



SHIPMENT DETAILS AND SERVICES





Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA



Railway Shipment



Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354

+97 150 696 1828 \ +97 155 676 1828 dubai@ecephan.com



 $(44" \times 44")$ 

Àsia



# We are good innovators

We do it **fast** 

&

We do it accurately





# **EHSP 800 LEVEL-1 ISOLATION GOWN**

AAMI Level-1 Isolation Gown (PP+PE or alternative SMS Material)





PP+PE (35 GSM) Exceeding Level - I

OR

35-40 GSM SMS Fabric



ISO 13485:2016 Certified



AAMI Level-I Tested





#### **GENERAL DESCRIPTION**

- Overlock Sewn
- Full Coverage with Overlapping
- The Neck Part Comes with an Overlapping Cloth, and the Back Neck Tie Comes with 2" and 4" Adjustable Velcro and has 1" Extended Overlapping Taps for Easy Gripping
- It is a Protective Procedure Cover Gown for Visitors or Standard Procedure Medical Units



AAMI Level-I Certified and Used for Basic Care with Minimal Risk



**Disposable** 



Sterility Non-Sterile



Gender Unisex



Fastening Type Adjustable Velcro Neck and Waist Ties



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL



**Packing Style** 50 PCs & 100 PCs

- Standard Procedure Medical Units
- Visitors

#### RECOMMENDED TO:

Low and Minimal Risk-Zones





Roadway Shipment

Airway Shipment



Railway Shipment

# e de la composición della composición della comp



SHIPMENT DETAILS AND SERVICES











CFR 2020





(1000 x 1200 mm) (40" x 48") Europe, North America

Standard -



1100 x 1100 mm) (44" x 44") Àsia



(800 x 1200 mm) (32" x 48") Europe





ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA



Swiss Office: ECEP HAN MEDIKAL TEK, SAN, TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



# **EHSP 810 ISOLATION GOWN**

AAMI Level-2 Isolation Gown (PP+PE or alternative SMS-SSMMS Material)





PP+PE (40 GSM) Exceeding Level - 2 Requirements

OR

SMS-SSMMS (40 GSM) Fabric



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms



AAMI Level-2 Tested





#### **GENERAL DESCRIPTION**

Over-lock Sewn

**EHSP-GOWNS** 

- Full Coverage with Overlapping
- The Neck Part Comes with an Overlapping Cloth, and the Back Neck Tie Comes with 2" and 4" Adjustable Velcro and Has 1" Extended Overlapping Taps for Easy Gripping
- Provides a Barrier Against Large Amount of Penetration through Splatter and Against Some Fluid Exposures through Soaking
- It is Used as a Protective Procedure Gown, Standard Isolation, and in Medical Units

#### RECOMMENDED TO:

- Low and Minimal Risk-Zones
- Standard Procedure Medical Units
- Visitors



AAMI Level-II Certified and Provides a Barrier Against Large Amount of Penetration



Disposable



Sterility Non-Sterile



Gender Unisex



Fastening Type Adjustable Velcro Neck and Waist Ties



**Private Labeling** Possible upon Request: Subject to MOQ.



S, M, L, XL, 2XL, 3XL



**Packing Style** 50 PCs & 100 PCs

# SHIPMENT DETAILS AND SERVICES



فشتت + اللك



Seaway Shipment



Roadway Shipment



Railway Shipment



EXW 2020











(1000 x 1200 mm) (40" x 48") Europe, North America

(1100 x 1100 mm)

Standard



 $(44" \times 44")$ Àsia



(800 x 1200 mm)



(32" x 48") Europe









ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101. SCOTTSDALE AZ 85260, ARIZONA USA














Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO, BOX NO. 39354





# **EHSP 820-3 SURGICAL GOWN**

AAMI Level-3 Surgical Gown (PP+PE or alternative SMS-SSMMS + Barrier)





#### **GENERAL DESCRIPTION**

- Over-lock Sewn
- Full Coverage with Overlapping
- The Neck Part Comes with an Overlapping Cloth, and the Back Neck Tie Comes with 2" and 4" Adjustable Velcro and has 1" Extended Overlapping Taps for Easy Gripping
- Used in Moderate Risk-Zones for Aerial Blood Draw and Inserting an Intravenous (IV) Line in Emergency Rooms, Medical Units, and Trauma Cases



AAMI Level-III Certified and Used in Moderate Risk-Zones



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender



Fastening Type Adjustable Velcro Neck and Waist Ties



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL



SHIPMENT DETAILS AND SERVICES

**Packing Style** 40 PCs & 80 PCs

- Emergency Rooms
- ICUs
- Trauma Cases









Seaway Shipment

Roadway Shipment

Railway Shipment



















(1000 x 1200 mm) (40" x 48") Europe, North America

(1100 x 1100 mm)

Standard



Àsia (800 x 1200 mm)

 $(44" \times 44")$ 



(32" x 48") Europe












AND



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354







# **EHSP 830-1 SURGICAL GOWN**

AAMI Level-3 Surgical Gown





Provides Barrier Against Some Fluid Exposures Through Soaking



AAMI Level-3 Tested



Provides a Barrier Against Large Amount of Penetration Through Splatter



40 SMS

25 GSM Barrier



#### GENERAL DESCRIPTION

- Full barrier Apron and Full Sleeves
- Over-lock Sewn
- Full Coverage with Overlapping
- The Neck Part Comes with an Overlapping Cloth, and the Back Neck Tie Comes with 2" and 4" Adjustable Velcro and has 1" Extended Overlapping Taps for Easy Gripping
- Used in Moderate Risk-Zones for Aerial Blood Draw and Inserting an Intravenous (IV) Line in Emergency Rooms, Medical Units, and Trauma Cases

AAMI Level-III Certified and Used in Moderate Risk Zones



**Disposable** 



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Adjustable Velcro Neck and Waist Ties



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL

+ 1 +



**Packing Style** 32 PCs & 65 PCs

- Healthcare Personnel During

#### **RECOMMENDED TO:**

Surgical Procedures

- ICUs

# Medium to High Risk-Zones

Roadway Shipment



Railway Shipment

Airway Shipment

Seaway Shipment

# G - 4 4 5

SHIPMENT DETAILS AND SERVICES











CFR 2020



(40" x 48") Europe, North America (1100 x 1100 mm)

Standard -



(44" x 44") Àsia



(800 x 1200 mm) (32" x 48")



Europe





ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK, SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND



swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354







# **EHSP 840-3 ISOLATION GOWN**

AAMI Level-2 Isolation Gown





40 GSM SMS 15 GSM Lamination 3 GSM PE (58 GSM)





Provides a Barrier Against Large Amount of Penetration Through Splatter



Provides Barrier Against Some Fluid Exposures Through Soaking



AAMI Level-2 Tested



#### **GENERAL DESCRIPTION**

- Full Barrier Apron Laminated
- Water and Fluid Splatter Protected with Full Sleeves
- Over-lock Sewn
- Full Coverage with Overlapping
- Neck Part Comes with Overlapping Cloth, and Back Neck Tie Comes with 2" and 4" Adjustable Velcro and has 1" Extended Overlapping Taps for Easy Gripping
- Provides a Barrier Against Large Amount of Penetration through Splatter and Against Some Fluid Exposures through Soaking. It is Used as a Protective Procedure Gown, Standard Isolation, and in Medical Units.

#### RECOMMENDED TO:

- Medium to High Risk-Zones
- Healthcare Personnel



AAMI Level-II Certified and Provides a Barrier Against Large Amount of Penetration



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Adjustable Velcro Neck and Waist Ties



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Free Size



**Packing Style** 50 PCs & 100 PCs

# SHIPMENT DETAILS AND SERVICES







TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA







Swiss Office: ECEP HAN MEDIKAL TEK, SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354







# **EHSP 850 CHEMOTHERAPY GOWN**

Disposable Chemotherapy Gown (ASTM 739-12), All Chemicals Tested





High Protection Against All Tested Chemicals



Barrier Against Large Amounts of Fluid Exposure over Long Periods



Barrier Against Blood Containing Virus, Pathogen & Infectious Diseases (Non-Airborne)



#### **GENERAL DESCRIPTION**

- Full Barrier Apron Laminated
- Water and Fluid Splatter Protected with Full sleeves
- Over-lock Sewn
- Full Coverage with Overlapping
- Neck Part Comes with Overlapping Cloth, and Back Neck Tie Comes with 2," and 4" Adjustable Velcro and has 1" Extended Overlapping Taps for Easy Gripping
- Microbial Barrier Protection to Protect the User Against Longtime Exposures
- Undergone ASTM F739-12 Standard Test Method for Permeation of Liquids and Gases Through Protective Clothing Materials Under Conditions of Continuous Contact
- Barrier Level Performance is Tested with Simulated Blood Containing a Virus, Pathogen Resistance, Infectious Diseases

#### **RECOMMENDED TO:**

- High Risk-Zones
- Healthcare Personnel



ASTM F739-12/ ASTM F1671 Certified and Used in High-Risk Situations

80 GSM

SMS WITH



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Adjustable Velcro Neck and Waist Ties



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL



**Packing Style** 40 PCs & 80 PCs

- Chemotherapy Purpose

## SHIPMENT DETAILS AND SERVICES









Railway Shipment











CFR 2020



(1000 x 1200 mm) (40" x 48") Europe, North America

1100 x 1100 mm)

Standard -



(44" x 44") Asia



(800 x 1200 mm) (32" x 48")

Europe



Innovation.Speed.Accuracy



TURKY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com

+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354

+97 150 696 1828 \ +97 155 676 1828 dubai@ecephan.com



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

+93 202 501 610 kabul@ecephan.com

#### **EHSP 860 L-1 & L-2 ISOLATION GOWN**

AAMI Level-1 and Level-2 Isolation Gown







Provides a Barrier Against Large Amount of Penetration Through Splatter



AAMI Level-2 Tested



AAMI Level-1 Tested



+10%



#### GENERAL DESCRIPTION

- Over-lock Sewn
- Full Coverage with Overlapping
- The Neck Part Comes with an Overlapping Cloth, and the Back Neck Tie Comes with 2" and 4" Adjustable Velcro and has 1" Extended Overlapping Taps for Easy Gripping
- Protective Procedure Cover Gown for Visitors or Standard Procedure Medical Units



AAMI Level-I & Level-II Certified & Used for Basic Care with Minimal Risk.



**Disposable** 



Sterility Non-Sterile



Gender Unisex



Fastening Type Adjustable Velcro Neck and Waist Ties



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL



**Packing Style** 50 PCs & 100 PCs

#### RECOMMENDED TO:

- Basic Care (Hospital Workers, Care Homes, or Low-Risk Wards)
- Hospital Visitors
- Hotels Cleanliness

## Airway Shipment Seaway Shipment











TURKEY OMCE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL











Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





+93 202 501 610 kabul@ecephan.com

#### **EHSP 900 L-4 REINFORCED SURGICAL GOWN**

AAMI Level-4 Reinforced Surgical Gown (PP+PE + Barrier + Tape Reinforcement or alternative SMS-SSMMS + Barrier + Tape Material)





PP+PE (55-60 GSM) 25 GSM BARRIER

Tape Reinforcement

OR

SMS-SSMMS (55-60 GSM)

25 GSM BARRIER

Tape Reinforcement



Provides Barrier Against Large Amount of Penetration Through Splatter



Barrier Against Blood Containing Virus, Pathogen & Infectious Diseases Using Phi-X174 Bacteriophage Penetration



Disposable

Gender

Unisex

AAMI Level-4 Tested





#### GENERAL DESCRIPTION

- Full Barrier Apron Laminated
- Water and Fluid Splatter Protected with Full sleeves
- Over-lock Sewn
- Full Coverage with Overlapping
- Neck Part Comes with Overlapping Cloth, Back Neck Tie Comes with 2," and 4" Adjustable Velcro and has 1" Extended Overlapping Taps for Easy Gripping
- Microbial Barrier Protection to Protect the User Against Longtime Contacting Blood Borne Pathogens
- ASTM F1671 Standard Test Method for the Resistance of Materials Used in Protective Clothing to Penetrate Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System

- Surgical Purpose



ANSI/AAMI Level -IV, ASTM F1671 Certified and Used in High-Risk Situations



Sterility Non-Sterile; Sterilized upon Request

S, M, L, XL, 2XL, 3XL

+ 11 + 11

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



**Fastening Type** 3" Knitted Cuff, Elasticated Back



Neck Velcro and Waist Tie-Back



**Private Labeling** Possible upon Request; Subject to MOQ.



SHIPMENT DETAILS AND SERVICES

€ - a i

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

**Packing Style** 32 PCs & 65 PCs

- High Risk-Zones
- Healthcare Personnel

Innovation.Speed.Accuracy

For Your Protection

#### **RECOMMENDED TO:**







ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA






Swiss Office:
ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND





DUBAI Office:
ECEP HAN INTERNATIONAL GENERAL TRADING LLC
TESLA BUSINESS CENTER OPAL TOWER,
SUIT NO. 90, 16th FLOOR 1610,
BUSINESS BAY PO. BOX NO. 39354





📮 Standard 🗅

(1000 x 1200 mm)

1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Europe

Àsia

ECEP HAN TRADING LLC

+93 202 501 610 kabul@ecephan.com



From blinking of the sunto hitting the hay we think about you & your HEALTH



#### EHSP 700 DISPOSABLE JACKET + TROUSER

Breathable, Lightweight & Flexible Material





PP+PE 55-60 GSM



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms





#### **GENERAL DESCRIPTION**

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases Micro-organisms such as Fluids, Aerosols, or Solid Dust Particles
- Not Stipulated by Seam Technology
- Tight Seam Sewn with all Seams Overlocked
- V-neck 1/4 of Full Sleeve
- Comfortable Wear for Relaxed Disposable Use at Low-Risk Environments



1/4 of Full Sleeve, Upper Shirt Pocket, and Two Trouser Pockets Plain Cuff without Elasticity



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type V-Neck, Disposable Waist-Tie Ban-Set



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL

+ 11 + .....



**Packing Style** 30 PCs & 60 PCs

- **Doctors**

- Patients

- Health Workers

#### **RECOMMENDED TO:**









Roadway Shipment



Railway Shipment

#### E - A 2



SHIPMENT DETAILS AND SERVICES











1100 x 1100 mm)  $(44" \times 44")$ Àsia

(40" x 48")

Standard S

(1000 x 1200 mm)

Europe, North America



(800 x 1200 mm) (32" x 48")

Europe





TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA



+1 480 656 5519 america@ecephan.com



Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



#### **EHSP 710 DOCTORS SCRUB**

SMS Hip-Lenght Scrub





40 GSM SMS Tested Fabric



Tested for Water Resistance-Impact Penetration



Tested for Water Resistance Hydrostatic Pressure





#### **GENERAL DESCRIPTION**

- SMS Disposable Doctors
- Breathable, Comfortable Sleeveless Pajamas
- Three Full-Size Pockets on the Shirt
- Elasticated Waistband
- Standard Color is Medical Blue, and any other Colors are Possible in Special Order
- AATCC 42:2017 Water Resistance-Impact Penetration Tested and AATC 127:2017 Tested



Breathable, Lightweight, and Flexible Material for Special Purposes Without Sleeve



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** V-Neck, Upper Shirt 3 Pocket & Elasticated Trouser Waist Band



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL

+ 11 +



SHIPMENT DETAILS AND SERVICES

**Packing Style** 40 PCs & 80 PCs

- Disposable for Nurses
- Disposable for Doctors
- Disposable for Healthcare Workers

#### **RECOMMENDED TO:**

### Airway Shipment



Roadway Shipment



Railway Shipment

#### **4** 4 2 INCOTERMS













(1000 x 1200 mm) (40" x 48") Europe, North America



Àsia (800 x 1200 mm)

(1100 x 1100 mm)

 $(44" \times 44")$ 



(32" x 48") Èurope



ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA



america@ecephan.com




Swiss Office:
ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



#### **EHSP 720 REUSABLE DOCTOR CLOTH**

Half Sleeve Scrub is made of 120 GSM Terry Cotton Fabric





Terry Cotton





#### **GENERAL DESCRIPTION**

- Durable
- Breathable
- Reusable
- Terry Cotton Fabric
- Comfortable
- Half Sleeve Scrub
- 1 Full-size Pocket on Shirt, Two Pockets on Pants
- Standard Color Medical Blue
- Other Colors are Possible in Special Orders



Hip-Length Reusable Doctor's Scrub Upper Shirt 3 Pocket, 2 Pants Pocket



Reusable



Sterility Non-Sterile



Gender Unisex



**Fastening Type** V-Neck, Waist Band for Tying



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL

+ 11 +



SHIPMENT DETAILS AND SERVICES

**Packing Style** 50 PCs & 100 PCs

#### **RECOMMENDED TO:**

- Reusable for Nurses
- Reusable for Doctors
- Reusable for Healthcare Workers

#### Airway Shipment Seaway Shipment



Roadway Shipment



Railway Shipment



EXW 2020











(1000 x 1200 mm) (40" x 48") Europe, North America

Standard 🗅



 $(44" \times 44")$ Àsia

(1100 x 1100 mm)



(800 x 1200 mm) (32" x 48") Èurope



ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA






Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



#### **EHSP 730 KNEE-LENGTH LAB COAT**

Light Blue 30 GSM Spunbond Light Weight 3 Pockets Knitted Lab Coat





30 GSM Spunbond Fabric



ISO 13485:2016 Certified





#### **GENERAL DESCRIPTION**

- Light Weight Fabric
- Light Blue
- White Knitted Cuffs
- Overlocked Seam
- Ultrasonic upon Request
- 3 Pockets
- V-Neck Tie Style
- Clamp on White Buttons



Full Sleeve with knitted cuff



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Clamp on Buttons Closure



**Private Labeling** Possible upon Request; Subject to MOQ.



S, M, L, XL, 2XL, 3XL

+ 11 +

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



SHIPMENT DETAILS AND SERVICES

**4** 4 2

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

**Packing Style** 50 PCs & 100 PCs

- Laboratory Workers
- Handling Biological, Chemical, or
- Handling Human Blood, Tissue, Cell Lines, or other Potentially Infectious Materials

- Radiological Materials and Waste

#### **RECOMMENDED TO:**



TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA












Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Standard -

(1000 x 1200 mm)

(1100 x 1100 mm)

(800 x 1200 mm) (32" x 48")

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

Àsia

Èurope

+93 202 501 610 kabul@ecephan.com



41

# Innovation.Speed.Accuracy

# HONESTY is our dogma we do best when there is no one to see



# ACCESSORIES



13" Tall Breathable Disposable Non-Slip Bootcover EN 14126 Certified Fabric





(PP+PE) 55-60 GSM Laminated Coverall



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms





#### **GENERAL DESCRIPTION**

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases, Micro-organisms such as Fluids, Aerosols, or Solid Dust Particles
- Not Stipulated by Seam Technology
- Multiple Thread Overlocked for Durable Usage
- EVA-Non Slip Soles makes it Super Comfortable
- Quality Product for Durable Use



Non-Slip EVA-Sole 13-15" Tall Bootcover, Elasticated on Cuffs, Full Ties around Cuff as well



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Elasticated over Cuffs & Tightening with Ankle Laces

+ 1 +

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



**Private Labeling** Possible upon Request; Subject to MOQ.



Size S, M, L



SHIPMENT DETAILS AND SERVICES

**4** 4 2

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

**Packing Style** 75 PCs & 150 PCs

- Mold Abatement
- Medical, Utilities, and Controlled Environments

#### **RECOMMENDED TO:**

- Pharmaceutical, Asbestos, Lead and

Innovation.Speed.Accuracy

 Food Processing, Manufacturing, Painting, Electronics, Automotive, and General Maintenance



TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkev@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA



Swiss Office:
ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

Standard S

(1000 x 1200 mm)

(1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Èurope

Àsia

+93 202 501 610 kabul@ecephan.com



Breathable, Lightweight, and Flexible Material. Used for Special Purposes and Special Environments, and Used in Low Risk Zones.





(PP+PE) 55-60 GSM



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms





#### GENERAL DESCRIPTION

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases, Micro-organisms such as Fluids, Aerosols, or Solid Dust **Particles**
- Not Stipulated by Seam Technology



4" Ankle Height Shoe-Cover Elasticated and Non-Slip Sole-less Plain Fabric Only.



**Disposable** 



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Elasticated over Ankles only



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Free Size



**Packing Style** 200 PCs & 400 PCs

#### **RECOMMENDED TO:**

- Healthcare Professionals to Minimize Contamination and Prevent Transmission of Micro-organisms between Staff and Patients
- Visitors of Hospitals and Clinics or Offices; to be not Used with Sharpie or Pointy Shoes

#### SHIPMENT DETAILS AND SERVICES





TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com


USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com

Swiss Office:
ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



Breathable, Lightweight, and Flexible Material Used for Special Purposes and Special Environments. Elasticated Overhead and Around the Face





(PP+PE) 55-60 GSM Laminated Coverall



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms





#### **GENERAL DESCRIPTION**

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases, Micro-organisms such as Fluids, Aerosols, or Solid Dust **Particles**
- Not Stipulated by Seam Technology



Covers Full Head, Back-Neck, Half Around Shoulders & Elasticated Around the Face



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Elasticated Around the Face



Private Labeling Possible upon Request; Subject to MOQ.



Standard Free Size

**4** + **1** + **1** 

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



**Packing Style** 200 PCs & 400 PCs

- Healthcare Professionals to Minimize Contamination and Prevent Transmission of Micro-organisms
- Medical Surgical Gowns
- Responders

#### **RECOMMENDED TO:**

- between Staff and Patients Used for Extra Protection under the
- Red-zone, Ambulances, First



Turkey Office: ECEP HAN MEDİKALTEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com




USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA



Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

🚍 Standard 🔼

(1000 x 1200 mm)

(1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Èurope

Àsia

+93 202 501 610 kabul@ecephan.com

46

SHIPMENT DETAILS AND SERVICES

🥌 🛶 🚜 gi

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020





#### GENERAL DESCRIPTION

- 20-25 Gram Comfortable Non-Woven Shoe-Cover
- Multiple Colors
- Elasticated for Extra Protection
- Packed in Dispenser Box of 100 PCs or 100 PCs/Bag
- Elasticated Over Feet, Easy to Wear, Super Light and Comfortable



Non-Woven Non-Sterilized Medical Graded 20-25 GSM Shoe-Cover



Disposable



Sterility Non-Sterile



Gender Unisex



**Fastening Type** Elasticated over Ankles Only



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Free Size



**Packing Style** 100 PCs / Plastic Bag

#### **RECOMMENDED TO:**

- Medical Hygienic Floors Protection Against Human Traffic
- Could be Used in Hospitals, Clinics or Offices
- Visitors of Hospitals and Clinics or Offices; to be not Used with Sharpie or Pointy Shoes

#### SHIPMENT DETAILS AND SERVICES

**4** 4 2 **4** + **1** + **1** Standard -INCOTERMS (1000 x 1200 mm) EXW 2020 Airway Shipment (40" x 48") FOB 2020 Europe, North America Seaway Shipment CIF 2020 (1100 x 1100 mm) (44" x 44") DAP 2020 Àsia Roadway Shipment DDP 2020 (800 x 1200 mm) (32" x 48") Railway Shipment CFR 2020 Èurope



Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA






Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

+93 202 501 610 kabul@ecephan.com





10-12 **GSM** 



Elasticated Over Head. Easy to Wear, Super Light and Comfortable





#### **GENERAL DESCRIPTION**

- 10-12 Gram Comfortable Bouffant
- Non-Woven Multiple Colors
- Elasticated to Protect your Hair
- Packed in Dispenser Box of 1000 PCs/Bulk Bag or 100 PCs/Bag
- Designed to Fit Adults with Long and Short Hair
- Elasticated Over Head, Easy to Wear, Super Light and Comfortable



Non-Woven, Medical Graded 10-12 GSM Bouffant



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Elasticated over Head



**Private Labeling** Possible upon Request; Subject to MOQ.



21" & 24" Longer Possible as well



Packing Style 100 PCs / Plastic Bag 1000 PCs / Bulk Bag

#### **RECOMMENDED TO:**

- Hygiene Protection in Medical Supplies, Industrial Safety
- Food-Processing Facilities, Kitchens, Beauty Salon and Houses

# Seaway Shipment

Airway Shipment

+ 1 +



Railway Shipment



SHIPMENT DETAILS AND SERVICES

EXW 2020

FOB 2020



DAP 2020

DDP 2020 CFR 2020



(40" x 48") Europe, North America



 $(44" \times 44")$ Àsia



(800 x 1200 mm) (32" x 48") Europe



ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com



Swiss Office:
ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI
37 ROUTE DE VIRELOUP
1239 COLLEX SWITZERLAND

+41 79 714 3173 swiss@ecephan.com




Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354

+97 150 696 1828 \ +97 155 676 1828 dubai@ecephan.com

Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

+93 202 501 610 kabul@ecephan.com

48



30-35 **GSM** 



Elasticated Over Head. Easy to Wear, Super Light and Comfortable





#### GENERAL DESCRIPTION

- 30-35 Gram Comfortable Bouffant
- Non-Woven Multiple Colors
- Elasticated to Protect your Hair
- Packed in Dispenser Box of 80 PCs or 100 PCs/Bag
- Designed to Fit Adults with Long and Short Hair
- 2" Wide Laces for Better Fastening
- Easy to Wear, Light and Comfortable



Non-Woven, Medical Graded 30-35 GSM Bouffant



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



Fastening Type Elasticated Around Head & Over Ear Laces



**Private Labeling** Possible upon Request; Subject to MOQ.



Standard Free Size

+ 1 +

Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment



SHIPMENT DETAILS AND SERVICES

**4** 4 2

INCOTERMS

EXW 2020

FOB 2020

CIF 2020

DAP 2020

DDP 2020

CFR 2020

Packing Style 300 PCs & 600 PCs / Dispenser Box 50 / Plastic Bag

#### **RECOMMENDED TO:**

Widely Used in Hospitals

Innovation.Speed.Accuracy

- Healthcare Professionals Extra Protections by Layering Against
- Particles and Low-Pressure Fluids

TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com




Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354



Standard -

(1000 x 1200 mm)

(1100 x 1100 mm)

(800 x 1200 mm)

Europe, North America

(40" x 48")

 $(44" \times 44")$ 

(32" x 48")

Èurope

Àsia

Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY









(PP+PE) 55-60 GSM Laminated Coverall



EN 14126 Standard Fabric



Protected Against Infectious Diseases & Micro-organisms





#### **GENERAL DESCRIPTION**

- EN 14126 Standard Defined Fabric
- Protected Against Infectious Diseases, Micro-organisms such as Fluids, Aerosols, or Solid Dust **Particles**
- Not Stipulated by Seam Technology
- Used in Medium-Risk Zones



Breathable, Lightweight, Flexible Material & Used for Special Purposes and Special Environments



Disposable



Sterility Non-Sterile; Sterilized upon Request



Gender Unisex



**Fastening Type** Elasticated at Wrist and Biceps



**Private Labeling** Possible upon Request; Subject to MOQ.



Size S, M, L



SHIPMENT DETAILS AND SERVICES

**Packing Style** 225 PCs & 450 PCs

#### **RECOMMENDED TO:**

- Healthcare Professionals for Extra Protection Against Pathogens and Patients Fluids by Possible Risks of Seam Punctures in the Medical Gowns
- Transferring Liquids and Solids in Medical Environment
- Industrial Usage

#### 🥌 🚤 🚜 💈

#### Airway Shipment

Seaway Shipment

+ 1 +



Roadway Shipment



Railway Shipment













(1000 x 1200 mm) (40" x 48") Europe, North America



(1100 x 1100 mm)  $(44" \times 44")$ Àsia

Standard -



(800 x 1200 mm) (32" x 48") Èurope





ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA






Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND



+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



#### **EHSP 2000 DISPOSABLE SPUNBOND BLANKET+PILLOW**

Non-Woven Blanket+Pillow Set







40 GSM Spunbond Blanket with 200 GSM Microfiber Filling + 40 GSM Spunbond Pillow with 550 GSM Microfiber



Sterility Non-Sterile; Sterilized upon Request



Packing Style Blanket+Pillow PE Pack Vacuumed



Gender Unisex



Size Blanket: 155x215 cm Pillow: 50x70 cm



Color Blue



**Private Labeling** Possible upon Request; Subject to MOQ.



Disposable

#### **GENERAL DESCRIPTION**

 Tensile Strength (GB/T 24218.3-2010): 62.2 (N), 41.0(N)

Elongation: 110.00%, 145.0(%)

Thickness: 0.37 (mm)

Water Vapor Transmission Rate: 7.67\*10<sup>3</sup>[g/(m<sup>2</sup>.24h)]

Water Repellency: 3-4 Class

Anti-electrostatic: Half life:3.8s

#### SHIPMENT DETAILS AND SERVICES

INCOTERMS

#### 🥌 👑 🚜 👩 **\*** + **1** + **....**











(1000 x 1200 mm) (40" x 48") Europe, North America

Standard



CIF 2020 DAP 2020



(1100 x 1100 mm)  $(44" \times 44")$ Àsia



Roadway Shipment

DDP 2020

CFR 2020



(800 x 1200 mm) (32" x 48") Èurope



TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND



+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





+93 202 501 610 kabul@ecephan.com

#### **DISPOSABLE FLEX TPE GLOVES**

Toxins free, tasteless and odorless, free of arcinogenic plasticizers as PVC





TPE material silicone quality elastomers a natural or synthetic rubber



New Generation Hybrid Technology



Not Allergic in Prolonged





#### **GENERAL DESCRIPTION**

- Perfect fit on hand and fingertips Very soft and skin-friendly Excellent tactile sensation
- Especially suitable for contact with greasy and oily foods Elastic and tear-resistant
- Odorless and tasteless
- %100 suitable for food contact Silicone and Latex-free
- Advantage of lower prices compared with alternative products Fast production process and short delivery times



TPE gloves are more smooth, more soft, more sensible, and more extensible.





**Disposable** 



**Private Labeling** Possible upon Request; Subject to MOQ.



SHIPMENT DETAILS AND SERVICES

**Packing Style** 50 PCs & 100 PCs

Non-Sterile

Sterility



**Fitting** Perfect Fitting Does not cause allergy.



Size S, M, L, XL

- Medical Areas, Hair & Beauty
- **Applications**

#### **RECOMMENDED TO:**

- Cleaning, Housekeeping
- Gastronomy
- Restaurants, Hotels
- Industrial Areas, Agricultural Areas



Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment





















CFR 2020



(1000 x 1200 mm) (40" x 48") Europe, North America

🚍 Standard 🔼



(1100 x 1100 mm) (44" x 44") Asia



(800 x 1200 mm) (32" x 48")



Europe





TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com

+41 79 714 3173 swiss@ecephan.com

37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND

Swiss Office



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354 ECEP HAN MEDIKAL TEK, SAN, TIC LTD. STI

+97 150 696 1828 \ +97 155 676 1828 dubai@ecephan.com



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



Toxins free, tasteless and odorless, free of arcinogenic plasticizers as PVC





TPE material silicone quality elastomers a natural or synthetic rubber



New Generation Hybrid Technology



Not Allergic in Prolonged



Powder-Free TPE Gloves





#### **GENERAL DESCRIPTION**

- Perfect fit on hand and fingertips Very soft and skin-friendly Excellent tactile sensation
- Especially suitable for contact with greasy and oily foods Elastic and tear-resistant
- Odorless and tasteless
- %100 suitable for food contact Silicone and Latex-free
- Advantage of lower prices compared with alternative products Fast production process and short delivery times



TPE gloves are more smooth, more soft, more sensible, and more extensible.



Disposable



Sterility Non-Sterile



Gender



**Fitting** Perfect Fitting Does not cause allergy.



**Private Labeling** Possible upon Request; Subject to MOQ.



Size S, M, L, XL



SHIPMENT DETAILS AND SERVICES

**4 4 3** 

**Packing Style** 50 PCs & 100 PCs

#### **RECOMMENDED TO:**

- Gastronomy
- Restaurants, Hotels
- Medical Areas, Hair & Beauty

- Cleaning, Housekeeping
- **Applications**
- Industrial Areas, Agricultural Areas



Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment





































(1000 x 1200 mm) (40" x 48") Europe, North America

Standard -



(1100 x 1100 mm) (44" x 44") Asia



(800 x 1200 mm) (32" x 48")

Europe





TURKEY O'MCE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com



+41 79 714 3173 swiss@ecephan.com

37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND

Swiss Office



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY







TPE material silicone quality elastomers a natural or synthetic rubber

> Worlinger ELDIVEN

> > Woods ELDIVEN

MAGO ELDIVEN



New Generation Hybrid Technology



Not Allergic in Prolonged



Powder-Free TPE Gloves





#### **GENERAL DESCRIPTION**

- Perfect fit on hand and fingertips Very soft and skin-friendly Excellent tactile sensation
- Especially suitable for contact with greasy and oily foods Elastic and tear-resistant
- Odorless and tasteless
- %100 suitable for food contact Silicone and Latex-free
- Advantage of lower prices compared with alternative products Fast production process and short delivery times



TPE gloves are more smooth, more soft, more sensible, and more extensible.



Wadie

ELDIVEN

Disposable



Sterility Non-Sterile



Gender



**Fitting** Perfect Fitting Does not cause allergy.



**Private Labeling** Possible upon Request; Subject to MOQ.



Size S, M, L, XL



SHIPMENT DETAILS AND SERVICES

**4 4 3** 

**Packing Style** 50 PCs & 100 PCs

#### **RECOMMENDED TO:**

- Gastronomy
- Restaurants, Hotels
- Medical Areas, Hair & Beauty
- **Applications**
- Industrial Areas, Agricultural Areas

- Cleaning, Housekeeping



Airway Shipment

Seaway Shipment

Roadway Shipment

Railway Shipment























CFR 2020



(1000 x 1200 mm) (40" x 48") Europe, North America

Standard -



(1100 x 1100 mm) (44" x 44") Asia



(800 x 1200 mm) (32" x 48")







+90 212 852 3447 turkey@ecephan.com



ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA







ECEP HAN MEDIKAL TEK, SAN, TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND

+41 79 714 3173 swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



#### **EHSP 2200 LEVEL 2 PATIENT KIT FULL PACKAGE**

Ecep Han
Trading for good.

The Kit Includes Five Products

#### **PACKAGE CONTENTS:**

- EHSP-560 Bouffant (1 Pc)
- EHSP-870 Level 2 Gown (1 Pc)
- EHSP-300 N95/FFP2 Mask (1 Pc)
- EHSP-550 Disposable Slipper (1 Pair)
- EHSP-540 Shoecover (1 Pair)




#### EHSP-560 10 GSM Bouffant 21" White 1 Pc

- 10-12 GSM Comfortable Bouffant
- Elasticated to Protect your Hair
- · Designed to Fit Adults with Long and Short Hair
- Elasticated Over Head, Easy to Wear, Super Light and Comfortable



#### EHSP-300 N95/FFP2 Mask 1 Pc

- NIOSH Tested N95 Respirator
- High-Efficiency Filtration > 95%<99.97%
- EN 149:2001 + A1:2010 Tested
- Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059 E13274-7:2019[%] Flow Rate 95 I/min Tested for Penetration of Paraffin Oil Mist



- Adjustable Velcro Neck and Waist Ties 1" Extended Overlapping Taps for Easy Gripping
- EN 14126 Standard Fabric
- Overlock and Ultrasonic Sewn



#### EHSP-540 SS 10 GSM Shoecover 1 Pair

- 20-25 Gram Comfortable
- Non-Woven Shoe-Cover
- Elasticated for Extra Protection
- Easy to Wear, Super Light and Comfortable



#### EHSP-550 Disposable Slipper 1 Pair

- UV Sterilized 50-60 GSM SB Fabric
- 50-60 gr Laminated SB innersole
- 3.5 mm Sole
- Size (40-47 fits all).





Turkey Office: ECEP HAN MEDİKALTEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA








37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND












#### EHSP 2300 LEVEL 1 PATIENT KIT BASIC PACKAGE

The Kit Includes Two Products



#### **PACKAGE CONTENTS:**

- EHSP-860 Level 1 Gown (1 Pc)
- EHSP-550 Disposable Slipper (1 Pair)





#### EHSP-860 Level 1 Gown 1 Pc

- 25-35 GSM SS Fabric
- MD (32) CD(45-75) +10%
- AAMI Level-1 & 2 Tested
- Adjustable Velcro Neck and Waist Ties 1" Extended Overlapping Taps for Easy Gripping
- Overlock or Ultrasonic Sewn
- Elasticated Overlock Wrist



#### EHSP-550 Disposable Slipper 1 Pair

- UV Sterilized 50-60 GSM SB Fabric
- 50-60 gr Laminated SB innersole
- 3.5 mm Sole
- Size (40-47 fits all).





TURKY OTTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL



turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA



america@ecephan.com



Swiss Office: ECEP HAN MEDIKAL TEK, SAN, TIC LTD, STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





DUDAI OMCE: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



#### **PACKAGE CONTENTS:**

- EHSP-860 Level 1 Gown (1 Pc)
- EHSP-540 Shoecover (1 Pair)
- EHSP-300 N95/FFP2 Mask (1 Pc)



#### EHSP-300 N95/FFP2 Mask 1 Pc

- NIOSH Tested N95 Respirator
- High-Efficiency Filtration > 95%<99.97%</li>
- EN 149:2001 + A1:2010 Tested
- Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059 E13274-7:2019[%] Flow Rate 95 I/min Tested for Penetration of Paraffin Oil Mist

#### EHSP-860 Level 1 Gown 1 Pc

- 25-35 GSM SS Fabric
- MD (32) CD(45-75) +10%
- AAMI Level-1 & 2 Tested
- Adjustable Velcro Neck and Waist Ties 1" Extended Overlapping Taps for Easy Gripping
- Overlock or Ultrasonic Sewn
- **Elasticated Overlock Wrist**



#### EHSP-540 SS 10 GSM Shoecover 1 Pair

- 20-25 Gram Comfortable
- Non-Woven Shoe-Cover
- Elasticated for Extra Protection
- Easy to Wear, Super Light & Comfortable





TURKEY O'MCE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkev@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA







Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND

ECEP HAN INTERNATIONAL GENERAL TRADING LLC



Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

+93 202 501 610 kabul@ecephan.com

# COME ON COME ALL OF WORKING WITH

ECEP HAN **ALL FOR** GOOD



# HYGIENE PRODUCTS



#### EHSP 1030, 1050 HAND SANITIZER WIPES









The Biocidal Product Type-I Hand Sanitizer Wipes are Used for Hand and Skin Disinfection. It Does not Need Water and with 78.5% Ethanol Kills Germs and Disinfects your Hands.

#### **INGREDIENTS**



#### **Aqua Water**

CAS Number:7732-18-5

EINECS/ELICS Number: 231-791-12

Function: Solvent

Concentration (%) (w/w): Up to 100



#### Ethanol

CAS Number: 64-17-5

EINECS/ELICS Number: 200-578-6 Function: Active Ingredient/Antibacterial Agent

SEA Classification: Flammable; H225 Eye Damage 2

H319, Spc Target Organ Tox

Concentration (%) (w/w): 78,5



#### Glycerin

CAS Number: 56-81-5

EINECS/ELICS Number: 200-289-5

Function: Humectant

Concentration (%) (w/w): 1,5



#### Aloe Barbadensis Leaf Juice

CAS Number: 85507-69-3/94349-62-9

EINECS/ELICS Number: 287-390-8/305-181-2

Function: Humectant

Concentration (%) (w/w): 0,1

**Colorless Transparent Liquid** 



4 ml Each Wipe



Odor



Characteristic pH Value @20°: 6,00-8,00



**Private Labeling** Possible upon Request; Subject to MOQ.



**Easy Tearable Sachet** 



**Packing Style** 24 PCs, and Individually Packed Bottles

#### SHIPMENT DETAILS AND SERVICES







+90 212 852 3447

TURKEY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA






Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





+93 202 501 610



#### **EHSP 1060** HAND SANITIZER STICKS









The Biocidal Product Type-I Hand Sanitizer is Used for Hand and Skin Disinfection. It Does not Need Water and with 78.5% Ethanol Kills Germs and Disinfects your Hands.

#### **INGREDIENTS**



#### **Agua Water**

CAS Number:7732-18-5

EINECS/ELICS Number: 231-791-12

Function: Solvent

Concentration (%) (w/w): Up to 100



#### Ethanol

CAS Number: 64-17-5

EINECS/ELICS Number: 200-578-6

Function: Active Ingredient/Antibacterial Agent SEA Classification: Flammable; H225 Eye Damage 2

H319, Spc Target Organ Tox

Concentration (%) (w/w): 78,5



#### Glycerin

CAS Number: 56-81-5

EINECS/ELICS Number: 200-289-5

Function: Humectant

Concentration (%) (w/w): 1,5



#### Aloe Barbadensis Leaf Juice

CAS Number: 85507-69-3/94349-62-9

EINECS/ELICS Number: 287-390-8/305-181-2

Function: Humectant

Concentration (%) (w/w): 0,1



**Colorless Transparent Liquid** 



4 mL Each Stick



#### Odour

Characteristic

pH Value @20°: 6,00-8,00



Private Labeling Possible upon Request; Subject to MOQ.



**Easy Tearable Sachet** 



**Packing Style** 10 PCs, Individually Packed in a Box

#### SHIPMENT DETAILS AND SERVICES

DDP 2020

CFR 2020





Railway Shipment



(800 x 1200 mm) (32" x 48")

(44" x 44")



Europe



TOIREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447 turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





+93 202 501 610



#### EHSP 1130, 1160, 1170 HAND SANITIZER BOTTLE









The Biocidal Product Type-I Hand Sanitizer is Used for Hand and Skin Disinfection. It Does not Need Water and with 78.5% Ethanol Kills Germs and Disinfects your Hands.

#### **INGREDIENTS**



#### **Aqua Water**

CAS Number:7732-18-5

EINECS/ELICS Number: 231-791-12

Function: Solvent

Concentration (%) (w/w): Up to 100



#### Ethanol

CAS Number: 64-17-5

EINECS/ELICS Number: 200-578-6

Function: Active Ingredient/Antibacterial Agent SEA Classification: Flammable; H225 Eye Damage 2

H319, Spc Target Organ Tox

Concentration (%) (w/w): 78,5



#### Glycerin

CAS Number: 56-81-5

EINECS/ELICS Number: 200-289-5

Function: Humectant

Concentration (%) (w/w): 1,5



#### Aloe Barbadensis Leaf Juice

CAS Number: 85507-69-3/94349-62-9

EINECS/ELICS Number: 287-390-8/305-181-2

Function: Humectant

Concentration (%) (w/w): 0,1



#### **Colorless Transparent Liquid**



50 ml, 500 ml, 1000 ml



#### Odour

Characteristic

pH Value @20°: 6,00-8,00



#### Private Labeling Possible upon Request;

Subject to MOQ.



**Packing Style** 24 PCs, and Individually Packed Bottles



Flip Top and Spray Type





#### +

EXW 2020

FOB 2020

CIF 2020

#### INCOTERMS

SHIPMENT DETAILS AND SERVICES



#### Standard (1000 x 1200 mm)

Airway Shipment





Seaway Shipment



Roadway Shipment

Railway Shipment



DDP 2020 CFR 2020



1100 x 1100 mm) (44" x 44")

Europe, North America

(40" x 48")



(800 x 1200 mm) (32" x 48")

Europe







TURKY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447

turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com







Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



#### **EHSP 1180** HAND SANITIZER REFILL GALLON







The Biocidal Product Type-I Hand Sanitizer is Used for Hand and Skin Disinfection. It Does not Need Water and with 78.5% Ethanol Kills Germs and Disinfects your Hands.

#### **INGREDIENTS**



#### **Aqua Water**

CAS Number:7732-18-5

EINECS/ELICS Number: 231-791-12

Function: Solvent

Concentration (%) (w/w): Up to 100



#### Ethanol

CAS Number: 64-17-5

EINECS/ELICS Number: 200-578-6

Function: Active Ingredient/Antibacterial Agent SEA Classification: Flammable; H225 Eye Damage 2

H319, Spc Target Organ Tox

Concentration (%) (w/w): 78,5



#### Glycerin

CAS Number: 56-81-5

EINECS/ELICS Number: 200-289-5

Function: Humectant

Concentration (%) (w/w): 1,5



#### Aloe Barbadensis Leaf Juice

CAS Number: 85507-69-3/94349-62-9

EINECS/ELICS Number: 287-390-8/305-181-2

Function: Humectant

Concentration (%) (w/w): 0,1



Odour

For Refill

Characteristic

**Colorless Transparent Liquid** 

pH Value @20°: 6,00-8,00



5000 mL



Private Labeling Possible upon Request; Subject to MOQ.



**Packing Style** 5 Liters (1.32 US Liquid Gallons) Single Packing















SHIPMENT DETAILS AND SERVICES





Airway Shipment





Seaway Shipment





Roadway Shipment Railway Shipment

EXW 2020 FOB 2020





DDP 2020 CFR 2020

(1000 x 1200 mm) (40" x 48") Europe, North America



(1100 x 1100 mm) (44" x 44")



(800 x 1200 mm) (32" x 48") Europe





TUREY OTICE: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL

+90 212 852 3447

turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA






Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND



swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



#### **EHSP 1190 HAND SANITIZER 20 LITERS**







The Biocidal Product Type-I Hand Sanitizer is Used for Hand and Skin Disinfection. It Does not Need Water and with 78.5% Ethanol Kills Germs and Disinfects your Hands.

#### **INGREDIENTS**



#### **Aqua Water**

CAS Number:7732-18-5 EINECS/ELICS Number: 231-791-12

Function: Solvent

Concentration (%) (w/w): Up to 100



#### Ethanol

CAS Number: 64-17-5

EINECS/ELICS Number: 200-578-6 Function: Active Ingredient/Antibacterial Agent

SEA Classification: Flammable; H225 Eye Damage 2

H319, Spc Target Organ Tox

Concentration (%) (w/w): 78,5



#### Glycerin

CAS Number: 56-81-5

EINECS/ELICS Number: 200-289-5

Function: Humectant

Concentration (%) (w/w): 1,5



#### Aloe Barbadensis Leaf Juice

CAS Number: 85507-69-3/94349-62-9

EINECS/ELICS Number: 287-390-8/305-181-2

Function: Humectant

Concentration (%) (w/w): 0,1

ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL



turkey@ecephan.com



USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA

+1 480 656 5519 america@ecephan.com




Swiss Office: ECEP HAN MEDIKAL TEK, SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND



swiss@ecephan.com



Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY



kabul@ecephan.com



**Colorless Transparent Liquid** 

#### Odour

Characteristic

pH Value @20°: 6,00-8,00



Gallon Openning



20 Liters



**Private Labeling** Possible upon Request; Subject to MOQ.




Packing Style Individually Packed

#### SHIPMENT DETAILS AND SERVICES







EHSP-GUIDE GOWN

















#### AAMI GOWN BARRIER LEVELS OF EHSP

The Association for the Advancement of Medical Instrumentation (AAMI) standards are designed to help medical professionals for the safe use of medical devices. AAMI introduced the voluntary standard ANSI/AAMI PB70:2012, Liquid Barrier Performance and Classification of Protective Apparel and Drapes Intended for Use in Health Care Facilities, to determine key identification measures for the appropriate selection of protective apparel and drapes, such as medical gowns, for use in healthcare facilities. The AAMI gown classification stems from four levels of barrier performance.



Minimal Fluid Barrier Protection: Used for daily patient care, when there is little to no risk of fluid exposure. Generally not used in the OR



**Moderate Fluid Barrier Protection:** Used for the widest range of surgical procedures, where the risk of fluid exposure is moderate.



Minimal to Low Fluid Barrier Protection: Used when there's only a slight risk of fluid exposure, for minimally invasive surgical procedures such as lumps and bumps.



**Highest Fluid and Microbial Barrier Protection:** Provides protection against blood-borne pathogens in critical zones; used for long, fluid-intensive procedures.

#### STANDARDS FOR ANSI/AAMI PB70:2012 GOWNS LEVEL 1 LEVEL 2 LEVEL 3 LEVEL 4 Moderate **BARRIER PROTECTION** Minimal Low High Liquid and Viral **RESISTANCE MEASURES** Liquid Penetration Liquid Penetration **Liquid Penetration** Penetration AATCC42 AATCC42 **TEST ASTM F1671** AATCC42 AATCC127 AATCC127 Spray Impact Spray Impact **TEST CRITERIA** ≤ 1.0g ≤ 1.0g Water Impact **Pass** ≤ 4.5g Hydrostatic Pressure Hvdrostatic Pressure ≥ 50cm ≥ 20cm Spray Impact 4% Spray Impact 4% ACCEPTABLE QUALITY LEVEL 4% 4% Hydrostatic Pr. 4% Hydrostatic Pr. 4%

The requirements of the European Union's EN13795 standard for the barrier capacity of surgical gowns mainly include resistance to dry and wet microbes and resistance to liquid penetration. Surgical gowns are divided into standard surgical gowns and reinforced surgical gowns according to different purposes. Ordinary surgical gowns are used for routine surgery with less bleeding, and reinforced surgical clothes are used for surgery with longer operation time and more bleeding, such as orthopedic surgery and gastrointestinal surgery.

#### STANDARDS FOR EN 13795 SURGICAL CLOTHES

QUALITY INDEX	STANDARD SURGICAL GOWN		REINFORCED SURGICAL GOWN	
	KEY AREA	SUB-KEY AREA	KEY AREA	SUB-KEY AREA
RESISTANCE TO MICROBIAL PENETRATION (DRY STATE) (LOG10CFU)	Not Required	≤2	Not Required	≤2
RESISTANCE TO MICROBIAL PENETRATION (WET) (BI)	≥2.8	Not Required	6.0	Not Required
RESISTANCE TO LIQUID PENETRATION (CMH2O)	≥20	≥10	≥100	≥10

#### **KNOW YOUR GOWNS**

#### **DEFINITION OF KEY TERMS:**



Resistance to microbial penetration (dry state) (log10CFU) represents the bacterial liquid with a concentration of 108 CFU/g, and shake the bacteria for 30 min. For the purpose of this standard,  $log10CFU \le 2$  means a maximum of 300 CFU.




**Resistance to microbial penetration (wet)** (BI) represents the lowest significant difference (LSD) of BI reached 0.98 in the 95% confidence interval, which is the lowest value between the two materials, and the difference is not within 0.98BI. BI=6.0, Indicates no penetration in this standard.



**Key area** Strictly maintain the sterile area during operation, including the front part of the waist, the upper part, the arms.



Sub-key area Areas with less sterility during operation, including the back, front and lower waist.





#### **EHSP 800 Level-1 Isolation Gown**

- PP+PE (35 GSM) Exceeding Level I
- AAMI Level-1 Tested
- Adjustable Velcro Neck and Waist Ties
- 1" Extended Overlapping Taps for Easy Gripping
- Overlock Sewn
- Neck Part Comes with an Overlapping Cloth, & the Back Neck Tie Comes with 2" & 4" Adjustable Velcro



#### **EHSP 820 Level-3 Isolation Gown**

- PP+PE (50-60 GSM)
- AAMI Level-3 Tested
- Adjustable Velcro Neck and Waist Ties
- 1" Extended Overlapping Taps for Easy Gripping
- Overlock Sewn
- Neck Part Comes with an Overlapping Cloth, & the Back Neck Tie Comes with 2" & 4" Adjustable Velcro



#### **EHSP 810 Level-2 Isolation Gown**

- PP+PE (40 GSM) Exceeding Level-2 requirements
- AAMI Level-2 Tested
- Adjustable Velcro Neck and Waist Ties
- 1" Extended Overlapping Taps for Easy Gripping
- Overlock Sewn EN 14126 Standard Fabric
- Neck Part Comes with an Overlapping Cloth, & the Back Neck Tie Comes with 2" & 4" Adjustable Velcro



#### **EHSP 830 Level-3 Surgical Isolation Gown**

- 40 SMS + 25 GSM Barrier (65 GSM)
- AAMI Level-3 Tested
- Adjustable Velcro Neck and Waist Ties
- 1" Extended Overlapping Taps for Easy Gripping
- Overlock Sewn
- Neck Part Comes with an Overlapping Cloth, & the Back Neck Tie Comes with 2" & 4" Adjustable Velcro



#### EHSP 840 Level-2 Isolation Gown

- 40 GSM SMS + 15 GSM Lamination + 3 GSM PE (58 GSM)
- AAMI Level-2 Tested
- Full Barrier Apron Laminated & Water & Fluid Splatter Protected with Full Sleeves
- Adjustable Velcro Neck and Waist Ties
- 1" Extended Overlapping Taps for Easy Gripping
- Overlock Sewn
- Neck Part Comes with an Overlapping Cloth, & the
- Back Neck Tie Comes with 2" & 4" Adjustable Velcro



#### EHSP 860 Level-1 & 2 Isolation Gown

- MD (32) CD(45-75) +10%
- AAMI Level-1 & 2 Tested
- Adjustable Velcro Neck and Waist Ties
- 1" Extended Overlapping Taps for Easy Gripping
- Overlock Sewn
- Neck Part Comes with an Overlapping Cloth, & the Back Neck Tie Comes with 2" & 4" Adjustable Velcro



#### **EHSP 850 Chemotherapy Gown**

- 80 GSM SMS with Lamination
- ASTM 739-12 & ASTM F1671 Certified All Chemicals Tested
- Full Barrier Apron Laminated & Water & Fluid Splatter Protected with Full Sleeves
- Adjustable Velcro Neck and Waist Ties
- 1" Extended Overlapping Taps for Easy Gripping
- Overlock Sewn
- Microbial Barrier Protection to Protect the Wearer Against Longtime Exposures




#### EHSP 900 Level-4 Reinforced Surgical Gown

- PP+PE (55-60 GSM) + 25 GSM BARRIER
   + Tape Reinforcement
- AAMI Level-4 Tested
- Full Barrier Apron Laminated & Water & Fluid Splatter Protected with Full Sleeves
- 3" Knitted Cuff, Elasticated Back Neck Velcro and Waist Tie-Back
- Overlock Sewn Microbial Barrier Protection to Protect the Wearer Against Longtime Exposures

#### **GOWN SIZE GUIDE**

GOWN PARTS	Small	Medium	Large	X-Large	2X-Large	3X-Large	
Gown Length	110/43	112/44	114/45	127/50	129/51	131/52	
Chests Width	120/47	130/51	135/53	140/55	145/57	150/59	Level 1
Sleeve Length	52.5/20.5	55/21.5	57/22.5	63/24.8	64/25	65/25.5	
Gown Length	110/43	112/44	114/45	135/53	137/54	139/55	
Chests Width	105/41	115/45	125/49	127/50	129/51	140/55	Level 2 & 3
Sleeve Length	52.5/20.5	55/21.5	57/22.5	63/24.8	64/25	65/25.5	
Gown Length	110/43	112/44	114/45	127/50	129/51	131/52	
Chests Width	125/49	135/53	140/55	150/59	155/61	160/63	Level 4 & Chemo Gown
Sleeve Length	52.5/20.5	55/21.5	57/22.5	63/24.8	64/25	65/25.5	down

Unit: cm/inch



70

#### PACKAGING INFORMATION



**Note:** This packaging information is mostly economical according to our estimates, however we are flexible with customers demands on different packaging requirements.

Description	Quantity	Non-Sterile	Sterile	Description	Quantity	Non-Sterile	Sterile
EHSP 800 Level				EHSP 830 Leve			
S	Idv Pck	EHSP800-NS-S	EHSP800-ST-S	S	ldv Pck	EHSP830-NS-S	EHSP830-ST-S
M	Idv Pck	EHSP800-NS-M	EHSP800-ST-M	M	ldv Pck	EHSP830-NS-M	EHSP830-ST-M
L	Idv Pck	EHSP800-NS-L	EHSP800-ST-L	L	ldv Pck	EHSP830-NS-L	EHSP830-ST-L
XL	Idv Pck	EHSP800-NS-X	EHSP800-ST-X	XL	ldv Pck	EHSP830-NS-X	EHSP830-ST-X
2XL	Idv Pck	EHSP800-NS-2X	EHSP800-ST-2X	2XL	ldv Pck	EHSP830-NS-2X	EHSP830-ST-2X
3XL	ldv Pck	EHSP800-NS-3X	EHSP800-ST-3X	3XL	ldv Pck	EHSP830-NS-3X	EHSP830-ST-3X
EHSP 810 Level	-2 Isolation			EHSP 850 Cher	notherapy G		
S	Idv Pck	EHSP810-NS-S	EHSP810-ST-S	S	ldv Pck	EHSP850-NS-S	EHSP850-ST-S
M	Idv Pck	EHSP810-NS-M	EHSP810-ST-M	M	ldv Pck	EHSP850-NS-M	EHSP850-ST-M
L	Idv Pck	EHSP810-NS-L	EHSP810-ST-L	L	ldv Pck	EHSP850-NS-L	EHSP850-ST-L
XL	Idv Pck	EHSP810-NS-X	EHSP810-ST-X	XL	ldv Pck	EHSP850-NS-X	EHSP850-ST-X
2XL	Idv Pck	EHSP810-NS-2X	EHSP810-ST-2X	2XL	ldv Pck	EHSP850-NS-2X	EHSP850-ST-2X
3XL	ldv Pck	EHSP810-NS-3X	EHSP810-ST-3X	3XL	ldv Pck	EHSP850-NS-3X	EHSP850-ST-3X
EHSP 840 Level	-2 Isolation	Gown		EHSP 900 Leve	I-4 Reinforce	ed Surgical Gown (W	hite)
S	Idv Pck	EHSP840-NS-S	EHSP840-ST-S	S	ldv Pck	EHSP900-NS-WE-S	EHSP900-ST-WE-S
M	Idv Pck	EHSP840-NS-M	EHSP840-ST-M	M	ldv Pck	EHSP900-NS-WE-M	EHSP900-ST-WE-N
L	Idv Pck	EHSP840-NS-L	EHSP840-ST-L	L	ldv Pck	EHSP900-NS-WE-L	EHSP900-ST-WE-L
XL	Idv Pck	EHSP840-NS-X	EHSP840-ST-X	XL	ldv Pck	EHSP900-NS-WE-X	EHSP900-ST-WE-X
2XL	Idv Pck	EHSP840-NS-2X	EHSP840-ST-2X	2XL	ldv Pck	EHSP900-NS-WE-2X	EHSP900-ST-WE-2
3XL	ldv Pck	EHSP840-NS-3X	EHSP840-ST-3X	3XL	ldv Pck	EHSP900-NS-WE-3X	EHSP900-ST-WE-3
EHSP 860 Level	-1 & 2 Isolat			EHSP 900 Leve	I-4 Reinforce	ed Surgical Gown (Bl	ue)
S	Idv Pck	EHSP860-NS-S	EHSP860-ST-S	S	ldv Pck	EHSP900-NS-BE-S	EHSP900-ST-BE-S
M	Idv Pck	EHSP860-NS-M	EHSP860-ST-M	M	ldv Pck	EHSP900-NS-BE-M	EHSP900-ST-BE-M
L	Idv Pck	EHSP860-NS-L	EHSP860-ST-L	L	ldv Pck	EHSP900-NS-BE-L	EHSP900-ST-BE-L
XL	Idv Pck	EHSP860-NS-X	EHSP860-ST-X	XL	ldv Pck	EHSP900-NS-BE-X	EHSP900-ST-BE-X
2XL	Idv Pck	EHSP860-NS-2X	EHSP860-ST-2X	2XL	ldv Pck	EHSP900-NS-BE-2X	EHSP900-ST-BE-2>
3XL	ldv Pck	EHSP860-NS-3X	EHSP860-ST-3X	3XL	ldv Pck	EHSP900-NS-BE-3X	EHSP900-ST-BE-3>
EHSP 820 Level	-3 Isolation						
S	ldv Pck	EHSP820-NS-WE-S				EXPERIENC	E OUR
		ELIODOGO NIO MENA	ELIOPOGO OT ME M				



#### To Order Please Contact

Idv Pck

ldv Pck

Idv Pck

Idv Pck

Idv Pck

EHSP 820 Level-3 Isolation Gown (Blue)

Idv Pck

Idv Pck

Idv Pck

Idv Pck

Idv Pck

Idv Pck



M

XL

2XL

3XL

S

M

XL

2XL

3XL

contracting@ecephan.com



+90 543 822 1717



www.ecephan.com



USA Office ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA





ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND





Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354



LEVEL 1

LEVEL 2

LEVEL 3

LEVEL 4

Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

**AAMI BARRIER LEVELS** 





Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLETI SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL



EHSP820-NS-WE-M EHSP820-ST-WE-M

EHSP820-NS-WE-X EHSP820-ST-WE-X

EHSP820-NS-WE-2X EHSP820-ST-WE-2X

EHSP820-NS-WE-3X EHSP820-ST-WE-3X

EHSP820-NS-BE-M EHSP820-ST-BE-M

EHSP820-NS-BE-2X EHSP820-ST-BE-2X

EHSP820-NS-BE-3X EHSP820-ST-BE-3X

EHSP820-ST-WE-L

EHSP820-ST-BE-S

EHSP820-ST-BE-L

EHSP820-ST-BE-X

EHSP820-NS-WE-L

EHSP820-NS-BE-S

EHSP820-NS-BE-L

EHSP820-NS-BE-X





















## EHSP-GUIDE M

#### ASTM FACE MASK PERFORMANCE LEVELS OF EHSP

Choosing a suitable mask in its proper place is one of the critical points for your PPE protocols. The American Society for Testing and Materials Standards (ASTM) is an internationally recognized organization that develops voluntary guidance for various industries and products to aid in quality control, product safety, and recommended usage and application.

ASTM standard F2100-11 provides specifications for face masks. ASTM levels for the face masks are classified into three levels of protection (low, moderate, high) to help clinicians decide which mask is appropriate for the task at hand. The European equivalent of ASTM Levels are classified in Types; Level 1 is equivalent to Type I, Level 2 is equivalent to Type II, Level 3 is equivalent to Type IIR, and N95 is equivalent to FFP2.

LVL1

Level 1 masks are suitable for brief

examinations, exposing radio-graphs, and cleaning tasks.



Level 2 masks are preferable for

procedures that involve a moderate level of aerosols, such as hand instrumentation and sealants.



Level 3 masks have a high level of protection and are

used for procedures involving high levels of aerosols such as ultrasonic scaling, surgical procedures, and crown preparation.

Medical face mask materials covered under ASTM specifications are designated as one or more of the following performance levels based on the barrier performance properties of the materials used only.

#### STANDARDS FOR ASTM F2100-11 LEVELS & EN 14683:2019 TYPES

	LEVEL 1	TYPE I	LEVEL 2	TYPE II	LEVEL 3	TYPE IIR
FLUID RESISTANCE	80 mmHg	Not Required	120 mmHg	Not Required	160 mmHg	@≥ 16 kPa > 120 mmHg
BFE	≥ 95%	≥ 95%	≥ 98%	≥ 98%	≥ 98%	≥ 98%
PFE, @ 0.1 MICRON	≥ 95%	Not Required	≥ 98%	Not Required	≥ 98%	Not Required
DELTA P	< 5.0 mm H <sub>2</sub> O/cm²	< 40.0 Pa/cm²	< 6.0 mm H <sub>2</sub> O/cm <sup>2</sup>	< 40.0 Pa/cm²	< 6.0 mm H <sub>2</sub> O/cm <sup>2</sup>	< 60.0 Pa/cm²
FLAME SPREAD	Class 1	2007/47/EC MDD 93/42/EEC	Class 1	2007/47/EC MDD 93/42/EEC	Class 1	2007/47/EC MDD 93/42/EEC

#### **DEFINITION OF KEY TERMS:**



**FLUID RESISTANCE** represents the mask's resistance to penetration by synthetic blood under pressure (mmHg). It measures the ability of a mask's material construction to minimize fluids from traveling through the material and potentially coming into contact with the wearer. The higher the fluid resistance (filtration), the better the protection.



**BFE** (Bacterial Filtration Efficiency) represents the percentage of aerosol particulates filtered at a size of 3 microns. It is the measure of the efficiency of the mask in filtering bacteria passing through it.



MASK SELECTION MADE EASY

**PFE** (Sub-micron Particulate Filtration Efficiency) represents the percentage of sub-micron particulates filtered at 0.1 microns. PFE is the measure of the efficiency of the mask in filtering particles passing through it. The size of the particles filtered is critical.




**DELTA P** (Differential Pressure) represents the pressure drop across the mask or resistance to air flow in mmH<sub>2</sub>O/cm<sup>2</sup> (European equivalent is shown as Pa/cm<sup>2</sup>). This determines breathing resistance – the higher the Delta P, the less the breathability, but the better the filtration.



**FLAME SPREAD** is a ranking derived by laboratory standard test methodology of a material's propensity to burn rapidly and spread flames. According to ASTM standards it is divided into classes and the European equivalent is EC directives.

#### **KNOW YOUR MASKS**



- The extra-long nose-bridge, enclosed aluminum nose piece forms a strong seal to increase protection.
- Seamless crimping makes it comfortable to wear.
- Meltblown filter layer for protection against viruses, bacteria and airborne particles.
- The elastic neckband provides ease in wearing and comfort in removal from the face during no-risk zones.
- Hydrophobic spunbond outer layer for maximum protection against fluids.
- 8 Anti splash and anti splatter fog free eye shield.
- The elasticated earloop is durable, breathable and comfortable to apply.
- 6 Soft and absorbent spunbond inner layer.
- 9 Vapor barrier to prevent fogging.





Masks Feature Guide  EHSP S-Mask, Earloop Style, Headloop Style, Tie-On, and procedure masks.  Compliant with ASTM and EN Standards.	Aluminum Adjustable Nose Bridge	Elastic Neckband	Elastic Earloop	Seamless Crimping	Spunbond Outer Layer	Spunbond Inner Layer	Meltblown Filter Layer	Fog Free Eye Shield	Vapor Barrier
EHSP 410 Level 1 (Type I) Disposable S-Mask	<b>/</b>	<b>/</b>	<b>/</b>	<b>/</b>	<b>/</b>	<b>/</b>	<b>/</b>		<b>/</b>
EHSP 410 Level 1 (Type I) Disposable Face Mask			<b>/</b>	/	/	<b>/</b>	<b>/</b>		
EHSP 415 Level 2 (Type II) Disposable S-Mask - Kiddo (Children)	✓×	/	/	<b>/</b>	/	/	/		
EHSP 415 Level 2 (Type II) Disposable Mask - Kiddo (Children)	✓×		/	<b>/</b>	/	/	/		<b>/</b>
EHSP 420 Level 2 (Type II) Disposable S-Mask	/	<b>/</b>	/	<b>/</b>	/	/	/		<b>/</b>
EHSP 420 Level 2 (Type II) Disposable Face Mask			/	/	/	/	/		
EHSP 430 Level 3 (Type IIR) Disposable S-Mask	/	/	/	/	/	/	/		/
EHSP 430 Level 3 (Type IIR) Disposable Face Mask	/		/	<b>/</b>	/	/	/		<b>/</b>
EHSP 440 Level 3 (Type IIR) Surgical Tie-On Mask	/			<b>/</b>	<b>/</b>	/	<b>/</b>		<b>/</b>
EHSP 450 Level 3 (Type IIR) Procedure Face Mask with Eye Shield	/		/	/	/	/	/		
EHSP 300 N95 (FFP2) Respirator Face Mask (Earloop)	<b>/</b>		/		/		/		
EHSP 310 N95 (FFP2) Respirator Face Mask (Headloop Style)	<b>/</b>			<b>/</b>	/		/		
EHSP 320 N95 (FFP2) Respirator Face Mask (Earloop)	<b>/</b>		<b>/</b>	<b>/</b>	<b>/</b>		<b>/</b>		
EHSP 330 N95 (FFP2) Respirator Face Mask (Earloop)	<b>/</b>		<b>/</b>	<b>/</b>	<b>/</b>		<b>/</b>		

#### BENEFITS OF S-MASKS:



- While driving alone and at no-risk zones, we usually tend to remove our masks, and in the process, we forget in the car or lose our masks, but with the S-mask you don't need to worry about it.
- Our kids lose masks more often, but with EHSP S-mask for Kids, we have solved this problem. Your kids will not lose their masks during no-wear zones.
- We not only care about your safety but also about your comfort and our S-Masks will provide the highest level of comfort a Mask can provide for you and your cute kids.
- S-masks not only provide a high level of comfort, but it saves your money and the climate by saving no.s of masks.







## We are good innovators

## We do it fast

&

## We do it accurately



Prolonged use of face mask can sometimes cause discomfort to the wearer. The EHSP S-Mask gives the wearer the liberty to easily unwear the mask at no risk-zones without compulsion of putting it in a pocket or a purse. Our Face Mask range of products gives our customers to choose between S-Masks and normal masks when ordering

from FHSP Products.





#### EHSP 410 Level 1 (Type I) Disposable S-Mask

- BFE> 95%
- PFE> 95% @ 0.1micron
- Fluid Penetration Resistance (80 mmHg)
- Delta P < 5.0mm H<sub>2</sub>O/cm<sup>2</sup>
- Sterilized upon Request (EO-Strilized)
- UV-Disinfected
- 3PLY Rectangular S-Mask



#### EHSP 410 Level 1 (Type I) Disposable Face Mask

- BFE> 95%
- PFE> 95% @ 0.1micron
- Fluid Penetration Resistance (80 mmHg)
- Delta P < 5.0mm H<sub>2</sub>O/cm<sup>2</sup>
- Sterilized upon Request (EO-Strilized)
- UV-Disinfected
- 3PLY Rectangular Face Mask



#### EHSP 415 Level 2 (Type II) Disposable S-Mask - Kiddo (Children)

- BFE> 98%
- PFE> 98% @ 0.1micron
- Fluid Penetration Resistance (120 mmHg)
- Delta P < 6.0mm H<sub>2</sub>O/cm<sup>2</sup>
- Sterilized upon Request (EO-Strilized)
- UV-Disinfected
- 3PLY Rectangular S-Mask & (Normal Mask)



#### EHSP 420 Level 2 (Type II) Disposable S-Mask

- BFE> 98%
- PFE> 98% @ 0.1micron
- Fluid Penetration Resistance (120 mmHg)
- Delta P < 6.0mm H<sub>2</sub>O/cm<sup>2</sup>
- Sterilized upon Request (EO-Strilized)
- UV-Disinfected
- 3PLY Rectangular S-Mask



#### EHSP 420 Level 2 (Type II) **Disposable Face Mask**

- BFE> 98%
- PFE> 98% @ 0.1micron
- Fluid Penetration Resistance (120 mmHg)
- Delta P < 6.0mm H<sub>2</sub>O/cm<sup>2</sup>
- Sterilized upon Request (EO-Strilized)
- **UV-Disinfected**
- 3PLY Rectangular S-Mask



#### EHSP 430 Level 3 (Type IIR) Disposable S-Mask & (Normal Mask)

- BFE> 98%
- PFE> 98% @ 0.1micron
- Fluid Penetration Resistance (160 mmHg)
- Delta P < 6.0mm H<sub>2</sub>O/cm<sup>2</sup>
- Sterilized upon Request (EO-Strilized)
- UV-Disinfected
- 3PLY Rectangular S-Mask & (Normal Mask)



#### EHSP 440 Level 3 (Type IIR) Surgical Tie-On Mask

- BFE> 98%
- PFE> 98% @ 0.1micron
- Fluid Penetration Resistance (160 mmHg)
- Delta P < 6.0mm H<sub>2</sub>O/cm<sup>2</sup>
- Sterilized upon Request (EO-Strilized)
- **UV-Disinfected**
- 3PLY Face Mask with Tie-On headband



#### EHSP 450 Level 3 (Type IIR) Procedure Face Mask with Eye Shield

- BFE> 98%
- PFE> 98% @ 0.1micron
- Fluid Penetration Resistance (160 mmHg)
- Delta P < 6.0mm H<sub>2</sub>O/cm<sup>2</sup>
- Sterilized upon Request (EO-Strilized)
- **UV-Disinfected**
- 3PLY S-Mask with Anti-Fog Eye-Shield



#### EHSP 340 N95/FFP2 Washable Mask

- EN 14683+AC:2019 BFE 99%
- DELTA P =  $0.85 \text{ H}_2\text{O/cm}^2$
- ASTM E2149 Anti Microbiol Activity
- Reduction 99.99% (Silver Material Filter)
   Eco-Friendly
- Washable and Reusable FFP2 Valve



#### **EHSP-PPP05** Face Shield

- Anti-Fog Visor, Gripper Mechanism
- Polycarbonate Plus 1 Spare Visors
- Adjustable Rubber Band to Fit Head Size
- Prevents Splash, Spray & Transfer of Non-Corrosive Fluids
- Impact Resistance: 595 KJm²

#### N95 FACE MASKS

We offer N95 headloop style masks (usually used in hospitals and by healthcare professionals), N95 earloop style masks (usually used by public and healthcare professionals), and N95 eco-friendly washable and reusable N95 masks. We give you the option to choose different colors of your choice.





#### EHSP 300 N95 (FFP2) Respirator Face Mask

- NIOSH Tested N95 Respirator
- High-Efficiency Filtration > 95%<99.97%
- EN 149:2001 + A1:2010 Tested
- Tested for Filter Standards 42CFR:84
   Equivalent of TEB-APR-STP-0059
- E13274-7:2019[%] Flow Rate 95 I/min Tested for Penetration of Paraffin Oil Mist



#### EHSP 310 N95 (FFP2) Respirator Face Mask (Headloop Style)

- NIOSH Tested N95 Respirator
- High-Efficiency Filtration > 95%<99.97%</li>
- EN 149:2001 + A1:2010 Tested
- Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059
- E13274-7:2019[%] Flow Rate 95 I/min Tested for Penetration of Paraffin Oil Mist



#### EHSP 320 N95 (FFP2) Respirator Face Mask (Earloop)

- NIOSH Tested N95 Respirator
- High-Efficiency Filtration > 95%<99.97%
- EN 149:2001 + A1:2010 Tested
- Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059
- E13274-7:2019[%] Flow Rate 95 I/min Tested for Penetration of Paraffin Oil Mist



#### EHSP 330 N95 (FFP2) Respirator Face Mask (Earloop)

- NIOSH Tested N95 Respirator
- High-Efficiency Filtration > 95%<99.97%
- EN 149:2001 + A1:2010 Tested
- Tested for Filter Standards 42CFR:84 Equivalent of TEB-APR-STP-0059
- E13274-7:2019[%] Flow Rate 95 I/min Tested for Penetration of Paraffin Oil Mist

## EHSP-GUIDE MASK

#### PACKAGING INFORMATION



Note: This packaging information is mostly economical according to our estimates, however we are flexible with customers demands on different packaging requirements.

#### ORDERING INFORMATION

Description	Quantity	Ref.#
EHSP 300 N95 (FFP2) F	Respirator Face Mask (	Earloop)
White	50/Box	EHSP300WE
White	60/Box	EHSP300WE
White	ldv Pck	EHSP300WE
EHSP 310 N95 (FFP2	) Respirator Face Ma	ısk (Headloop)
Blue	50/Box	EHSP310BE
Pink	50/Box	EHSP310PK
White	50/Box	EHSP310WE
Dlook	EO/Day	FLICDATORIZ

EHSP 320 N95 (FFP2)	Respirator Face Ma	sk (Earloop)
Black	50/Box	EHSP310BK
White	50/Box	EHSP310WE
Pink	50/Box	EHSP310PK

EU25 350 M32 (LL55)	nespirator race ina	isk (⊑arioop)
Black	50/Box	EHSP320BK
Black	60/Box	EHSP320BK
Black	Idv Pck	EHSP320BK

EHSP 330 N95 (FFP2) Res	pirator Face M	ask (Earloop)
Blue	50/Box	EHSP330BE
Pink	50/Box	EHSP330PK
White	50/Box	EHSP330WE
Black	50/Box	EHSP330BK

EHSP 410 Level 1 (Type I) Disposable S-Mask					
Blue	50/Box	EHSP410-S-BE			
Pink	50/Box	EHSP410-S-PK			
White	50/Box	EHSP410-S-WE			
Black	50/Box	EHSP410-S-BK			
Purple	50/Box	EHSP410-S-PE			

EHSP 410 Level 1	(Type I) Disposable Fac	ce Mask
Blue	50/Box	EHSP410BE
Pink	50/Box	EHSP410PK
White	50/Box	EHSP410WE
Black	50/Box	EHSP410BK
Purple	50/Box	FHSP410PF

EHSP 415 Level 2 (Type II)	Disposable S-Ma	sk - Kiddo (Children)
Blue	50/Box	EHSP415-S-BE
Pink	50/Box	EHSP415-S-PK
Red	50/Box	EHSP415-S-RD
Purple	50/Box	EHSP415-S-PE

EHSP 415 Leve	I 2 (Type II)	Disposable Mask	- Kiddo	(Children)

Blue	50/Box	EHSP415BE
Pink	50/Box	EHSP415PK
Red	50/Box	EHSP415RD
Purple	50/Box	EHSP415PE

EHSP 420 Level 2 (Type II) Disposable S-Mask

**Description** 

	( · ) p · · · / - · · p · · · · · · · · · · · · · · · ·	
Blue	50/Box	EHSP420-S-BE
Pink	50/Box	EHSP420-S-PK
White	50/Box	EHSP420-S-WE
Black	50/Box	EHSP420-S-BK
Purple	50/Box	EHSP420-S-PE

Quantity

Ref.#

EHSP 420 Level 2 (Type II)	Disposable Fac	ce Mask
Blue	50/Box	EHSP420BE
Pink	50/Box	EHSP420PK
White	50/Box	EHSP420WE
Black	50/Box	EHSP420BK
Purple	50/Box	EHSP420PE

EHSP 430 Level 3 (Type IIR)	Disposable S	6-Mask
Blue	50/Box	EHSP430-S-BE
Green	50/Box	EHSP430-S-GN
White	50/Box	EHSP430-S-WE
Black	50/Box	EHSP430-S-BK
Purple	50/Box	EHSP430-S-PE

EHSP 430 Level 3 (Type	e IIR) Disposable	Face Mask
Blue	50/Box	EHSP430BE
Green	50/Box	EHSP430GN
White	50/Box	EHSP430WE
Black	50/Box	EHSP430BK
Purple	50/Box	EHSP430PE

EHSP 440 Level 3 (Type IIR	(a) Surgical Tie	-On Mask
Green	50/Box	EHSP440GN
White	50/Box	EHSP440WE
Blue	50/Box	EHSP440BE

EHSP 450 Level 3 (Type IIR) Prod	cedure Face N	lask with Eye Shield
Blue	50/Box	EHSP450BE
Green	50/Box	EHSP450GN
White	50/Box	EHSP450WE

50/Box

#### To Order Please Contact



sales@ecephan.com



+90 543 822 1717



www.ecephan.com





Black



LEVEL 1

LEVEL 2



Turkey Office: ECEP HAN MEDİKAL TEK. SAN. TİC LTD. ŞTİ ÇAKMAKLI MAH. DEVLET SOK NO:4/A, 34500 BÜYÜKÇEKMECE / İSTANBUL





USA Office: ECEP HAN NORTH AMERICA LLC 8260 E. GELDING DR. SUITE 101, SCOTTSDALE AZ 85260, ARIZONA USA







Swiss Office: ECEP HAN MEDIKAL TEK. SAN. TIC LTD. STI 37 ROUTE DE VIRELOUP 1239 COLLEX SWITZERLAND







Dubai Office: ECEP HAN INTERNATIONAL GENERAL TRADING LLC TESLA BUSINESS CENTER OPAL TOWER, SUIT NO. 90, 16th FLOOR 1610, BUSINESS BAY PO. BOX NO. 39354





Afghanistan: ECEP HAN TRADING LLC HOUSE NO 754, STREET NO. 14 DISTRICT 6, POSTAL CODE NO. 1006 KART-E-SAY

EHSP450BE



### CERTIFICATES AND TEST REPORTS



















## Sertifika Certificate

#### ECEP HAN MEDİKAL TEKSTİL SAN. VE TIC. LTD. STI.

AKÇABURGAZ MAHALLESİ AKÇABURGAZ CADDESİ NO 20 / 28 AKTİM 1 TİCARET VE İŞ MERKEZİ ESENYURTI İSTANBUL

DSR BELGELENDİRME tarafından denetlenmiş ve uygulamakta olduğu Kalite Yönetim Sisteminin

KAL ÜRÜNLERİN VE MEDİKAL EKİPMANLARIN KİMYASAL, BEŞERİ (İLAÇ, GIDA, MEDİKAL) VE HÜYENİK ÜRÜNLERİN ÜRETİMİ İTHALATI, İHRACATI VE ULUSLARARASI TİCARETİ

PRODUCTION, IMPORT, EXPORT AND INTERNAL TRADE OF MEDICAL PRODUCTS AND MEDICAL EQUIPMENT, OHEMICAL, HUMAN (PHARMACEUTICAL, FOOD, MEDICAL) AND HYGIENIC PRODUCTS

Sertifika No / Certificate No: QMS-21.04.338

05.08.2021 Sertifika Son Basım Tarihi Certificate Last Issue Date 3 Yıl/Years

04.08.2022 Sertifika Geçerlilik Tarihi Certificate Expiry Date

EA Kodu/EA 23-29









DSR TEKNÍK MUAYENE SERTÍFÍKASYON HÍZ. LTD. \$TÍ. Sk. Hells Beyaz Ofis 271 A Blok No.287 PENDÍK / ISTANBUL / TÜRKÍYE Tel: +98 216 506 06 35 ní: Infall/distribiglelendimus zan. 1r - Welt: www.distriblelendimus.com. 1r

## Sertifika Certificate

#### ECEP HAN MEDİKAL TEKSTİL SAN. VE TİC. LTD. STİ.

AKÇABURGAZ MAHALLESÎ AKÇABURGAZ CADDESÎ NO 20 / 28 AKTÎM 1 TÎCARET VE ÎŞ MERKEZÎ ESENYURTI ÎSTANBUL

DSR BELGELENDİRME tarafından denetlenmiş ve uygulamakta olduğu Çevre Yönetim Sisteminin

PRODUCTION, IMPORT, EXPORT AND INTERNAL TRADE OF MEDICAL PRODUCTS AND MEDICAL EQUIPMENT, OHEMICAL HUMAN (PHARMACEUTICAL FOOD, MEDICAL) AND HYGIENIC PRODUCTS

Sertifika No / Certificate No: EMS-21.04.338

03.05.2021

03.05.2021 3 Yil/Years

02.05.2022 Sertifika Geçerilik Tanhi Certificate Expiry Date EA Kodu/EA 23-29







OSR TEKNIK MUAYENE SERTIFIKASYON HIZ. LTD. \$TI.

Yeşilbağlar Man. Selvis Sa. Neisi Beyzo Olis 27 A bina No-220 PENDIK. ISTANBUL, YÜNCİYE Tel: +90 216 406 06 36
Email: info@dsrbajqelendirme,com.tr - Web: www.dsrbajqelendirme.com.tr

## QUALITY CERTIFICATIO



#### SERTIFIKA / CERTIFICATE

ECEP HAN MEDÎKAL TEKSTÎL SANAYÎ VE TÎCARET LİMİTED ŞİRKETİ

AKÇABURGAZ MAH. AKÇABURGAZ CAD. NO:20 AD ESENYURT / İSTANBUL / TURKEY

Kupulusunus <sup>44</sup>Tek kullanimlik cerrahi maske tek kullanimlik yöz maskesi tek kullanisilik cerrahi önlök, tek kullanimlik medikal önlök, tek kullanimlik eldiven tek kullanimlik tulum, tek kullanimlik miayene örtüleri, tek kullanimlik eldiven, tek kullanimlik tulum, tek kullanimlik bone, galos, sülkom maske, filtreli sülkön maske, koruyucu gözlök, tek kullanimlik sedye öriösö üretimi, satişi ve İHRACATI" Kapsamı için

For scope "DEPOSABLE SURGICAL MASK, DISPOSABLE FACE MASK, DISPOSABLE SURGICAL GOWNS, DISPOSABLE MEDICAL GOWNS, DISPOSABLE SURGICAL DRAPE, DISPOSABLE INSPECTION COVERS, DISPOSABLE GOWES, DISPOSABLE CONFERLLS, DISPOSABLE DISPOSABLE MASS. BLICON, MASK, FILTERED SILICON, MASK, PROTECTIVE GLASSES, DISPOSABLE STREET/FIER COVER PRODUCTION, SALES AND EXPORT

ISO 13485: 2016

#### TIBBİ CİHAZLAR KALİTE YÖNETİM SİSTEMİ Medical Devices Quality Management System


Kurduğunu ve uyguladığını belgelemekte ve NVA tarafından gerçekleştirilen denetim bu yönetim sisteminin yukarıda belirtilen standardın şartlanın karşıladığını doğrulamaktadır.

It certifies that it is established and implemented, and the audit performed by NVA it confirms that this management system meets the requirements of the following standard.

Sertifika Numarası / Certificate Number Sertifika Kodu / Certificate Code Sertifika Yayın Taribi / Certificate Issue Date Sertifika Gerefilik Taribi, Certificate Validity Date Sertifika Periyodu / Certificate Period

- : 20091805 : ECEP HAN MEDİKAL





retirin bir pakitite sünstüslistige ve götellim fatkilları zamanında yaşınlığır müddeteye bu halpe 1 yıl geçenlilir. NVA dissibil alteklinde geretil ilina ve yelkinin gödelmesene reğimen böyük hömellerde danil sorumluluk kabul dimeysceklir. Bu belbe VVA dirir ve kitindiğinde jakde edilmişinde



This Certificate has been awarded to

#### ECEP HAN MEDİKAL TEKSTİL SAN. VE TİC. LTD. ŞTİ AKÇABURGAZ MAH. AKÇABURGAZ CAD. NO.20/28 AKTİM 1 TİCARET VE İŞ MERKEZİ ESENYURT / İSTANBUL / TÜRKİYE

In recognition of the organization's Management System which complies with

#### ISO 45001:2018 (OHSMS)

The scope of activities covered by this certificate is defined below

PRODUCTION, IMPORT, EXPORT AND INTERNAL TRADE OF MEDICAL PRODUCTS AND MEDICAL EQUIPMENT, CHEMICAL, HUMAN (PHARMACEUTICAL, FOOD, MEDICAL) AND HYGIENIC PRODUCTS

MEDİKAL ÜRÜNLERİN VE M<mark>EDİKAL EKİP</mark>MANLARIN, KİMYASAL, BEŞERİ (İLAÇ, GIDA, MEDİKAL) VE HİJYENİK ÜRÜNLERİN ÜRETİMİ, İTHALATI, İHRACATI VE ULUSLARARASI TİCARETİ

Certificate Number: SISTURO05202126791

Date of Issue of Original Certificate: 03.05.2021 Date of Issue of latest certificate: 03.05.2021 Expiry Date: 02.05.2022



The Crganization's documentation and implementation has been reviewed an comply with the relevant standard rules. This certificate of Registration is bas evaluation of the mentioned scope given above. Organization is responsible for the responsibilities of the relevant standard rules. Any significant changes in the certification or standard referred above render this certificat This is an accredited certificate issued by SIS Certifications Pvt. Ltd. sanctioned fraternational Accreditation Services. 3000 Saturn Street State 100 Bres., California 92821

Corporate Office (SIS):- Plot No. 1539, 2nd Floor, Sector-4, Gurgaon -122001, Haryana Email us :- info@elscertifications.co.in, arun@siscertifications.us Call:- 0124-466 Web:- http://www.siscertifications.co.in, www.siscertifications The status of this certificate can be verified on "http://www.siscertifications.



IAF



#### ATTESTATION CERTIFICATE

Technical file of company mentioned below has been inspected and audit has been completed successfully Medical Devices 93/42/EEC have been taken as references for these processes

ECEP HAN MEDIKAL TEKSTIL SAN. VE TIC. LTD. STI Company Name

AKÇABURĞAZ MAH, AKÇABURĞAZ CAD, NO 20 / 28 AKTİM 1 TİCARET VE IŞ MERKEZİ ESENYURTI İSTANBUL MEDICAL DEVICES 93/42/EEC Company Address

Related Directives and Annex

Related Standards

Product Name

Product Brand/Model/Type

PROTECTIVE COVERALL! Coverall, Doctor S FINDIECTIVE COVERALL GOVERNI, Doctor Scrubs, Bouffant, Shoel Bootcover, Handover, Overales Produce, Nutrile Glovers, Overalls Produc Codes: EHSP-400, EHSP-4608, EHSP-4608, EHSP-4608, EHSP-4608, EHSP-4608, EHSP-4608, EHSP-4608, EHSP-4608, EHSP-4608, EHSP-4609, EHSP

Certificate Number Registration Date 05.08.2021

: CE-21.08.201

Expiry Date 04.08.2022







Technical file of company mentioned below has been inspected and audit has been completed successfully Medical Devices 93/42/EEC have been taken as neferences for these processes

ECEP HAN MEDIKAL TEKSTIL SAN, VE TIC, LTD, STI Company Name

: AKÇABURGAZ MAH. AKÇABURGAZ CAD. NO 20 / 28 AKTİM 1 TİCARET VE IŞ MERKEZİ ESENYURTI İSTANBUL. : MEDICAL DEVICES 93/42/EEC Company Address

Related Directives and Annex

Related Standards Product Name

: TS EN 13795-2 : EHSP : SURGICAL GOWN (Disposable Medic Cost Apron product codes: EHSP 800, EHSP-30-1, EHSP-460,

CE-21.08.201-1 Registration Date 05.08.2021

Expiry Date 04.08.2022

The CE mark si











#### CERTIFICATE OF REGISTRATION

The International Bureau of the World Intellectual Property Organization (WIPO) certifies that the indications appearing in the present certificate conform to the recording made in the International Register of Marks maintained under the Madrid Agreement and Protocol.

Reproduction of the mark in color under Rule 9.4(a)(vii)



1 613 894 Registration date April 14, 2021

Date next payment due April 14, 2031

Name and address of holder

Classification of 26.4; 27.5; 29.1 figurative elements List of goods and services NCL(11-2021)

ECEP HAN MEDIKAL TEKSTIL SANAYI VE TICARET LIMITED

Akcaburgaz Mah. Akcaburgaz Cad. 20 AD, Esenyurt İstanbul (Turkey)

3 Bleaching and cleaning preparations, detergents other than for use in manufacturing operations and for medical purposes, laundry bleach, fabric softeners for laundry use, stain removers, dishwasher detergents; perfumery; non-medicated cosmetics; fragrances; deodorants for personal use and animals; soaps; dental care preparations; dentifitices, denture polishes, tooth whitening preparations, mouth washes, not for medical purposes; abrasive preparations; emery cloth; sandpaper; pumice stone; abrasive pastes; polishing preparations for leather, vinyl, metal and wood, polishes and creams for leather, vinyl, metal and wood, polishes and creams for leather, vinyl, metal and wood, polishes and ereams for leather, vinyl, metal and wood, polishes and ereams for leather, vinyl, metal and wood, polishes and ereams for leather, vinyl, metal and wood, polishes and ereams for leather, vinyl, metal and wood, polishes and ereams for leather, vinyl, metal and veterinary preparations for medical purposes; chemical reagents for pharmaceutical and veterinary purposes; dietary supplements for pharmaceutical and veterinary purposes; dietary supplements; nutritional supplements; medical preparations for silmming purposes; food for babies; herbs and herbal beverages adapted for medicinal purposes; dental impression material, dental adhesives and material for repairing teeth; sanitary preparations for medical use; hygienic pads; hygienic tampons; plasters; materials for dressings; diapers made of paper and textiles for babies, adults and pets; preparations for destroying vermin; herbicides, fungicides, preparations for destroying vermin; herbicides, fungicides, preparations for destroying vermin; herbicides, fungicides, preparations for destroying vermin; herbicides, fungicides, preparations for destroying vermin; herbicides, fungicides.

WORLD STATE OF STATE

# **I QUALITY CERTIFICA**



#### SERTIFIKA / CERTIFICATE

ECEP HAN MEDİKAL TEKSTİL SANAYÎ VE TİCARET LIMITED SIRKETI

> AKCABURGAZ MAH. AKCABURGAZ CAD, NO:20 AD ESENYURT / ISTANBUL / TURKEY

Kuruluşunun "ALKOLSÜZ ANTISEPTIK MENDIL" Kapsamı için

For scope "NON-ALCOHOL ANTISEPTIC WIPES"

#### MS 1500:2009 HELAL / HALAL / حلال

Kurduğunu ve uyguladığını belgelemekte ve NVA tarafından gerçekleştirilen denetim bu yönetim sisteminin yukarıda belirtilen standardın şartlarını karşıladığını doğrulamaktadır.

It certifies that it is established and implemented, and the audit performed by NVA it confirms that this management system meets the requirements of the following sta

Sertifika Numarası / Certificate Number Sertifika Kodu i Certificate Code Sertifika Yayın Tarihi / Certificate Issue Date Sertifika Geçerlilik Tarihi / Certificate Valldity Date Sertifika Periyodu / Certificate Period : 20091601

ECEP HAN MEDİKAL 16.09.2020 : 16.09.2021 : 1 Vil / J Yo

kleri zamanında yapıldığı müddetpe bu belge 1 yıl geçerlidir. NVA denetim redmen büyük ilkimallarda dahil sonumbulur kabul atmayarasisir. Bu helsem















82007-2017 Regions Co







ingsir dan tegerar b





TEST / INSPECTION REPORT **EUROLAB LABORATORY SERVICES** 



2020121867

ECEP HAN MEDIKAL TEKS. SAN. LTD. STİ. / Akçaburgaz mah. Akçaburgaz cad.
No.20/28 Aktim 1 Ticaret ve iş merkezi Esenyurt-İSTANBUL
ECEP HAN NORTH AMERICA LLC/ 8260 E. Gelding Dr. Sulte 101, SCOTTSDALE Az

85260, Arizona- USA
37 ROUTE DE VIRELOUP/ 1239 Collex SWITZERLAND BUSSINES BAY/ Omniyat Building Tower PO. Box No:39354

Nikmal ABDULLAH 0212 852 3447 /+14806565519 / +41797143173 / +971567233175 **Contact Telephone** Contact e-mail: n.com / America Vecephen.com / Swiss Secephan.com /

Dubai@ecephan.com mple Accepted on: Report Date: 18.12.2020 Total number of pages:

EHSP-300 N95 Mask (Tekli Paket-50li Kutu)

EHSP-300 N95 Mask (60li Kutu)

TEST		METHOD	Specimen	RESULT
Respiratory protective devices - Filtering half masks to protect EN 149:2001	EHSP-300 N95 Mask (Tekli	PASS		
	EN 149:2001	Paket-50li Kutu)	FFP2	
	nst particles - Requirements, testing, marking	+A1:2009	EHSP-300 N95	PASS
testing, marking		Mask (60li Kutu)	FFP2	





Merkez Mh. Gençosman Cd, No 11 / A. GUNGÖREN / ISTANBUL Tal: 0212 702 20 10 Fax: 0212 909 21 10. Web: www.laboratuvar.com E-mail: info@laboratuvar.com



TEST / INSPECTION REPORT **EUROLAB LABORATORY SERVICES** 



2020121868

**Applicant**:

20201213696

ECEP HAN MEDÍKAL TEKS. SAN. LTD. ŞTİ. / Akçaburgat mah. Akçaburgat cad.
No.30/28 Aktım. 1 Ticaret ve iş merikezi Esanyuri. ISTANBUL

ECEP HAN NORTH AMERICA LLC/ 8260 E. Gelding Dr. Suite 101, SCOTTSDALE At

85260, Arizona- USA
37 ROUTE DE VIRELOUP/ 1239 Collex SWITZERLAND

BUSSINES BAY/ Omniyat Bullding Tower PD. Box No.39354 Nikmal ABDULLAH 0212 852 3447 / +14805565519 / +41797143173 / +971567233175

Contact Telephone: E-Mail: m.com / America@lecept Dubai@ecephan.com

Sample Accepted on: 03.12.2020 Report Date: Total number of pages 18 12 2020

Sample ID: EHSP-300 N95 Mask (Tekli Paket-50li Kutu)

EHSP-300 N95 Mask (60h Kutu)

١	TEST	METHOD	Specimen	RESULT
	Sodium Chloride (NaCl) Aerosol	NIOSH N95	EHSP-300 N95 Mask (Tekli Paket-50li Kutu)	PASS
	Test for 3D Filter Mask	NIOSH NSS	EHSP-300 N95 Mask (60li Kutu)	PASS



**EUROLAB** 

Report No:

Applicant:

Contact Person: Contact Telephone:

Sample Accepted on: Report Date: Total number of pages:

TEST

ce of Materials Used

Standard Specification for

in Medical Facial Masks

Determining the First Efficiency of Substances Used in Medical Facial Masks to Penetration Using Particles Using Latex Spheres

Sample ID:





5190243I802

2020121864

Meriez Mh. Gençosman Cd. No 11 / A. GUNGOREN / ISTANBLIL Tel: 0212 702 20 10 Fax: 0212 909 21 10 Web: www.laboraturar.com E-mail: Info@taboraturar.com

TEST/MUAYENE RAPORU

**EUROLAB LABORATUVAR HİZMETLERİ** 

ECEP HAN MEDIKAL TEKS. SAN. LTD. ŞTİ. / Akçaburgaz mah. Akçaburgaz cad. No:20/28 Aktim

VALUE

LIIVE Z

W98,45

1498,39

1 Ticaret ve is merkezi Esenyurt- ISTANBUL ECEP HAN NORTH AMERICA LLC/ 8260 E. Gelding Dr. Suife 101, SCOTTSDALE Az 85260,

37 ROUTE DE VIRELOUP/ 1239 Collex SWITZERLAND

METHOD SPECIMENS

@ecephan.com / Ame i@ecephan.com

ASTM F 2100

ASTM F 2299

03.12.2020

18.12.2020

BUSSINES BAY/ Omniyat Building Tower PG. Box No.39354 Nikmal ABDULLAH 0212 852 3447 /+14806565519 /+41797143173 /+971567233175

EHSP-430 Type IIR S-Mask (Yetişkin) / EHSP-430 Type IIR S-Mask (Cocuk)

EHSP-430 Type IIR S-Mask

(Yetişkin) EHSP-430 Typ IIR S-Mask

(Çocuk) EHSP-430 Type IIR S-Masi

> IIR 5-Mask (Çocuk)



TEST / INSPECTION REPORT **EUROLAB LABORATORY SERVICES** 



2020121863-2

ECEP HAN MEDIKAL TEKS. SAN, LTD. STI. / Akçaburgaz mah. Akçaburgaz cad. No.20/28 Aktim 1 Ticaret ve iş merkezi Esenyuri-ISTANBUL ECEP HAN NORTH AMERICA LLC/ 8260 E. Gelding Dr. Suite 101, SCOTTSDALE Az

37 ROUTE DE VIRELOUP/ 1239 Collex SWITZERLAND

BUSSINES BAY/ Omniyat Building Tower PO, 80x No:39354
Nikmal ABDULLAH
0212 852 3447 / +14806565519 / +41797143173 / +971567233175 **Contact Telephone** 

velkan@erohan.com / America@ecephan.com / Swiss@ecephan.com /
Dubai@ecephan.com

Contact e-mail:

Report Date: 18.12.2020 Total number of pages: 9 (Pg)

EHSP-415 TYPE IIR S-Mask (Çocuk)

TEST	МЕТНОВ	Specimen	RESULT
Medical face masks - Requirements and		EHSP-415 TYPE IIR S-Mask (Çocuk)	PASS
test methods	EN 14683+AC 2019		TYPE IIR











Hasan KUTLU



RESULTS

PASS

PASS

PASS

Hava Sanavdin

rikez Mh. Gengosman Cd. No 11 / A GONGOREN / ISTANBUL Tel: 0212 702 20 10 Fax: 0212 909 21 10 Wati www.laboratuvar.com E-mail: info@raboratuvar.com



#### TEST/MUAYENE RAPORU

#### **EUROLAB LABORATUVAR HİZMETLERİ**



Report No:

ECEP HAN MEDİKAL TEKS. SAN. LTD. ŞTİ. / Akçaburgaz mah. Akçaburgaz cad. No:20/28 Aktim Applicant:

1 Ticaret ve iş merkezi Esenyurt: ISTANBUL ECEP HAN NORTH AMERICA LLC/ 8260 E, Gelding Dr. Suitte 101, SCOTTSDALE Ar 85260.

37 ROUTE DE VIRELOUP! 1239 Colles SWITZERLAND BUSSINES BAY/ Omniyat Building Tower PQ: Box No:39354 Nikmal ABDULLAH

Contact Person :

0212 852 3447 / +14806565519 / +41797143173 / +971567233175 Contact Telephone:

Dubai@ecephan.com / Ame 03.12.2020 Sample Accepted on:

ort Date: 18.17.2020 Total number of pages: 7 (Pg)

Sample ID:

EHSP-430 Type IIR S-Mask (Yetişkin) / EHSP-430 Type IIR S-Mask (Çocuk)

	TEST	METHOD	SPECIMENS	VALUE	RESULTS	
	Standard Specification for	ASTM F 2100	EHSP-430 Type IIR S-Mask (Yetişkin)	Level 2	PASS	
	Performance of Materials Used in Medical Facial Masks	ASTM F 2100	EHSP-430 Type IIR S-Mask (Cocuk)	Level 2	PASS	
		Standard Test Method for Determining the First Efficiency of Substances Used in Medical	3	EHSP-430 Type IIR S-Mask (Yetişkin)	%98,45	PASS
	Facial Masks to Penetration Using Particles Using Latex Spheres	ASTM F 2299	EHSP-430 Type IIR S-Mask (Çocuk)	1698,39	PASS	



Hasan KUTLU



Hava Sariavdin

kez Mh. Gençosman Cd. No 11 / A. GÜNGÖREN / ISTANBUL Tel: 0212 702 20 10 Fax: 6212 909 21 10 Web: evew laboratuvar.com E-mail: info@laboratuvar.com



TEST / INSPECTION REPORT **EUROLAB LABORATORY SERVICES** 



2020121615

ECEP HAN MEDİKAL TEKS. SAN, LTD. ŞTİ. / Akçaburgaz mah. Akçaburgaz cad.
No:20/28 Aktim 1 Ticaret ve işi merkazi Esanyurt- İSTANBUL

ECEP HAN NORTH AMERICA LLC/ 8260 E. Gelding Dr. Suite 101, SCOTTSDALE Az

85260, Arizona-USA

37 ROUTE DE VIRELOUP/ 1239 Collex SWITZERLAND

BUSSINES BAY/ Omniyat Building Tower PO. Box No:39354 Nikmal ABDULLAH 0212 852 3447 / +14806565519 / +41797143173 / +971567233175 volkan@ecephan.com / America@ecephan.com / Swiss@ecephan.com / Dubai@ecephan.com E-Mail:

Sample Accepted on 20.11.2020

Report Date: 18.12.2020

Sample ID: **EHPS-850 DISPOSABLE ISOLATION GOWN** 

TEST	МЕТНОО	SPECIMEN	RESULT
Permeation Of Liquids And Gases Through Protective Clothing Materials Under Conditions Of Continuous Contact	ASTM F 739-12	EHPS-850 DISPOSABLE ISOLATION GOWN	PASS
Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System	ASTM F 1671	EHPS-850 DISPOSABLE ISOLATION GOWN	PASS







z Mh, Gengoeman C6, No 11 / A. GUNGOREN / ISTANBUL Tet, 0212 700 20 10 Fax: 0212 909 21 10 st; www.laboratunar.com E-mail: info@laboratunar.com



TEST / INSPECTION REPORT **EUROLAB LABORATORY SERVICES** 



2020121869
ECEP HAN MEDIKAL TEKS. SAN. LTD. \$TI. / Akçaburgaz mah. Akçaburgaz cad.
No.20/28 Aktim 1 Ticaret ve ili merkezi Esenyurt- ISTANBUL
ECEP HAN NORTH AMERICA LLC/ 8260 E. Geiding Dr. Suite 101, SCOTTSDALE Az

7 (Pg)

37 ROUTE DE VIRELOUP/ 1239 Collex SWITZERLAND

BUSSINES BAY/Omniyat Bullding Tower PO. Box No:39354 Nikmal ABDULLAH 0212 852 3447 / +14806565519 / +41797143173 / +971567233175 Contact Telepho Contact e-mail:

volkan@ecephan.com / America@ecephan.com / Swiss@ecep

Report Date: 18.12.2020

mple ID:

EHSP-900 Non Steril Surgical Reinfurced Isolation Gown

TEST	METHOD	Specimen	RESULT
Classification of barrier performance of surgical gowns, other protective apparel, surgical drapes and drape accessories	ANSI / AAMI PB 70:12	EHSP-900 Non Steril Surgical Reinfurced Isolation Gown	Level a
			PASS
Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System	ASTM F 1671	EHSP-900 Non Steril Surgical Reinfurced Isolation Gown	PASS







Merkez Mr. Gençosman Cd. No 11 / A GONGOREN / ISTANBUL. Tel: 0212 702 20 10 Fax: 0212 909 21 10 Web www.latocatives.com.



705 Front Street Toledo, OH 43605 hone: 419-693-5307 Fax: 419-691-0418

ECEP HAN MEDIKAL TEKSTIL VE TICARET LTD STI Akcaburgaz Mahallesi Akcaburgaz Caddesi NO:20 Suite 28 Aktim 1 Ticaret Ve Is Merkezi, Esenyurt/Istanbul Attn: Nick Hammond Rev 01 ECEP HAN NORTH AMERICA LLC 8260 E. Gelding Dr. Suite 101 Scottdale AZ, 85260 ARIZONA LISA

Revised Certificate

AAMI PB70 Liquid Barrier Performance and Classification

Test Article: EHSP 820-3 Gown

A total of thirty-two (32) specimens were tested from thirty-two (32) test articles. Specimens were chosen from the critical zones as described in AAMI PB70 for an isolation gown. Test specimens were subjected to the following tests:

AATCC 42 Water Resistance: Impact Penetration Test. AATCC 127 Water Resistance: Hydrostatic Pressure Test.

Based on the results of the testing as summarized in the attached reports, numbers 2004657 and 2004668, the product listed above was classified as AAMI PB70 Level 3.

Revision Statement: This certificate was revised to correct the Istanbul address and to add the USA address. The test results were not affected.

Record Storage: All raw data pertaining to this study will be maintained in the LexaMed archives for a minimum of 5 years.

Derison Duli





Page 1 of 1 FR100847-R0-052120











## Our Other Brands







www.bobosync.com

Lead the way









